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皮膚感作性試験代替法
Interleukin-8 Reporter Gene Assay
(IL-8 Luc assay)

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**JaCVAM statement on
the Interleukin-8 Reporter Gene Assay (IL-8 Luc assay)
Skin Sensitization Test Method**

At a meeting held on 25 June, 2020, at the National Institute of Health Sciences (NIHS) in Kanagawa, Japan, the Japanese Center for the Validation of Alternative Methods (JaCVAM) Regulatory Acceptance Board unanimously endorsed the following statement:

Proposal: Although it is possible to classify chemicals that yield positive results using this IL-8 Luc assay test method as skin sensitizers, it is difficult to classify accurately their skin sensitization strength. If a negative result is obtained, it is not possible to judge as a non-skin sensitizer. Therefore, the IL-8 Luc assay test method is not suitable for predicting skin sensitization potential on its own; in order to make a suitable assessment, the results of IL-8 Luc assay testing must be used with a thorough understanding of the properties of each test chemical in combination with other information as part of an integrated approach to testing and assessment (IATA). Furthermore, when using this test, thorough consideration must be given to the applicability domain and to the judgement of inconclusive chemicals.

This statement was prepared, following the review prepared by the Skin Sensitization Test JaCVAM Editorial Committee, to acknowledge that the results of a review and study by the JaCVAM Regulatory Acceptance Board have confirmed the usefulness of this assay.

Based on the above, we propose the IL-8 Luc assay test method as a useful means for estimating skin sensitization by regulatory agencies.


Yasuo Ohno
Chairperson

JaCVAM Regulatory Acceptance Board

Yoko Hirabayashi
Chairperson

JaCVAM Steering Committee

December 2, 2020

The JaCVAM Regulatory Acceptance Board was established by the JaCVAM Steering Committee, and is composed of nominees from the industry and academia.

This statement was endorsed by the following members of the JaCVAM Regulatory Acceptance Board:

Mr. Yasuo Ohno (Kihara Memorial Yokohama Foundation for the Advancement of Life Sciences): Chairperson

Ms. Yoko Hirabayashi (Center for Biological Safety and Research: CBSR, National Institute of Health Sciences: NIHS)

Mr. Morihiko Hirota (Japan Cosmetic Industry Association)

Mr. Yoshiaki Ikarashi (NIHS)

Mr. Takanori Ikeda (Japan Pharmaceutical Manufacturers Association)**

Mr. Noriyasu Imai (Japanese Society for Alternatives to Animal Experiments)

Mr. Kunifumi Inawaka (Japan Chemical Industry Association)

Mr. Tomoaki Inoue (Japanese Society of Immunotoxicology) *

Mr. Yuji Ishii (CBSR, NIHS)

Ms. Yumiko Iwase (Japan Pharmaceutical Manufacturers Association)*

Mr. Fumihiko Kubo (Pharmaceuticals and Medical Devices Agency)*

Mr. Kenichi Masumura (Japanese Environmental Mutagen Society)

Ms. Ruriko Nakamura (National Institute of Technology and Evaluation)

Mr. Akiyoshi Nishikawa (CBSR, NIHS/ Saiseikai Utsunomiya Hospital)

Mr. Jihei Nishimura (Pharmaceuticals and Medical Devices Agency)

Mr. Satoshi Numazawa (Japanese Society of Toxicology)

Ms. Keiko Yamamoto (Pharmaceuticals and Medical Devices Agency)**

Mr. Hiroo Yokozeki (Japanese Society for Cutaneous Immunology and Allergy)

Term: From 1st April 2018 to 30th September 2020

*: From 1st April 2018 to 31st March 2020

** : From 1st April 2020 to 30th September 2020

This statement was endorsed by the following members of the JaCVAM steering Committee after receiving the report from JaCVAM Regulatory Acceptance Board:

Ms. Yoko Hirabayashi (CBSR, NIHS): Chairperson
Mr. Osamu Fueki (Pharmaceuticals and Medical Devices Agency)
Mr. Yukihiro Goda (NIHS)
Mr. Akihiko Hirose (Division of Risk Assessment, CBSR, NIHS)
Mr. Koji Ishii (National Institute of Infectious Diseases)
Mr. Yasunari Kanda (Division of Pharmacology, CBSR, NIHS)
Mr. Satoshi Kitajima (Division of Toxicology, CBSR, NIHS)
Ms. Kumiko Ogawa (Division of Pathology, CBSR, NIHS)
Mr. Masahiro Sasaki (Ministry of Health, Labour and Welfare)
Mr. Keiichi Sugiyama (Division of Genetics and Mutagenesis, CBSR, NIHS)
Mr. Masahiro Takahata (Ministry of Health, Labour and Welfare)
Mr. Yuhji Taquahashi (Animal Management Section of the Division of Toxicology, CBSR, NIHS)
Mr. Masaaki Tsukano (Ministry of Health, Labour and Welfare)
Mr. Masahiko Yokota (Pharmaceuticals and Medical Devices Agency)
Mr. Takao Ashikaga (Division of Risk Assessment, CBSR, NIHS): Secretary
Mr. Hajime Kojima (Division of Risk Assessment, CBSR, NIHS): Secretary

皮膚感作性試験代替法
Interleukin-8 Reporter Gene Assay (IL-8 Luc assay)

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評価会議報告書	1
評価報告書	7
OECD KEY EVENT BASED TEST GUIDELINE 442E, IN VITRO SKIN SENSITISATION ASSAYS ADDRESSING THE KEY EVENT ON ACTIVATION OF DENDRITIC CELLS ON THE ADVERSE OUTCOME PATHWAY FOR SKIN SENSITISATION.....	37

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- 1) Takahashi T, Kimura Y, Saito R, Nakajima Y, Ohmiya Y, Yamasaki K, Aiba S. (2011) An in vitro test to screen skin sensitizers using a stable THP-1-derived IL-8 reporter cell line, THP-G8. *Toxicol Sci.* Dec;124:359-69.
- 2) OECD (2017). Validation report for the international validation study on the IL-8 Luc assay as a test evaluating the skin sensitizing potential of chemicals conducted by the IL-8 Luc Assay. Series on Testing and Assessment No. 267, ENV/JM/MONO (2017)19. Organization for Economic Cooperation and Development, Paris. Available at: <https://www.oecd.org/env/ehs/testing/series-testing-assessment-publications-number.htm>
- 3) Report of the peer review panel for the IL-8 Luciferase (IL-8 Luc) assay for in vitro skin sensitization. Series on testing and assessment. No.268 ENV/JM/MONO(2017)20. Organization for Economic Cooperation and Development, Paris. Available at: <http://www.oecd.org/env/ehs/testing/series-testing-assessment-publications-number.htm>
- 4) OECD TG442E (2018) In vitro skin sensitization assays addressing the key event on activation of dendritic cells on the adverse outcome pathway for skin sensitisation. Accessible at: http://www.oecd-ilibrary.org/environment/oecd-guidelines-for-the-testing-of-chemicals-section-4-health-effects_20745788
- 5) JaCVAM 皮膚感作性資料編纂委員会：皮膚感作性試験報告書 Interleukin-8 Reporter Gene Assay (IL-8 Luc assay) (2019年11月1日)
- 6) Kimura Y, Watanabe M, Suzuki N, Iwaki T, Yamakage K, Saito K, Nakajima Y, Fujimura C, Ohmiya Y, Omori T, Kojima H, Aiba S. (2018) The performance of an in vitro skin sensitisation test, IL-8 Luc assay (OECD442E), and the integrated approach with direct peptide reactive assay (DPRA), *J. Toxicol. Sci.*, 43 (12), 741-749.
- 7) Weber FC1, Németh T, Csepregi JZ, Dudeck A, Roers A, Ozsvári B, Oswald E, Puskás LG, Jakob T, Mócsai A, Martin SF. (2015) Neutrophils are required for both the sensitization and elicitation phase of contact hypersensitivity, *J. Exp. Med.*, 212 (1), 15-22.
- 8) Liu T, Zhang L, Joo D, Sun SC. (2017) NF- κ B signaling in inflammation, *Signal Transduct Target Ther.* 2, 17023.

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AOP:	Adverse Outcome Pathway
AP-1:	Activator Protein 1
ARE:	Antioxidant/electrophile Response Element
ARE-Nrf2:	ARE-Nuclear factor-erythroid 2-related factor 2
BrdU:	5-Bromo-2-deoxyuridine
CRO:	Contract Research Organization
DC:	Dendritic Cell
DMSO:	Dimethyl sulfoxide
DPRA:	Direct Peptide Reactivity Assay
EU:	European Union
FBS:	Fetal Bovine Serum
GAPDH:	GlycerAldehyde 3-Phosphate Dehydrogenase
IATA:	Integrated Approaches to Testing and Assessment
h-CLAT:	Human Cell Line Activation test
ICATM:	International Cooperation on Alternative Test Methods
IL-8:	Interleukin-8
IL-8 Luc Assay:	Interleukin-8 Reporter Gene Assay
JaCVAM:	Japanese Centre for the Validation of Alternatives Methods
KE:	Key Event
KeratinoSens TM :	The ARE-Nrf2 luciferase KeratinoSens TM test method
LA:	Lactic acid
LLNA:	Local Lymph Node Assay
LuSens:	The ARE-Nrf2 luciferase LuSens test method
MTA:	Material Transfer Agreement
4-NBB:	4-NitroBenzyl Bromide
NFkB:	Nuclear Factor kB
OECD:	Organisation for Economic Co-operation and Development
QSAR:	Quantitative Structure-Activity Relationship
REACH:	Registration, Evaluation, Authorization and Restriction of Chemicals
SLO:	Stable Luciferase Orange
SLR:	Stable Luciferase Red
TG:	Test Guideline
U-SENS TM :	U937 cell line activation Test
UN GHS:	United Nations Globally Harmonized System of Classification and Labelling of Chemicals
VMT:	Validation Management Team

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1st run	2nd run	3rd run	4th run	Final prediction
Positive	Positive	-	-	Positive
	Negative	Positive	-	Positive
		Negative	Positive	Positive
		Negative	Negative	Supposed negative
Negative	Positive	Positive	-	Positive
		Negative	Positive	Positive
		Negative	Negative	Supposed negative
	Negative	Positive	Positive	Positive
		Negative	Negative	Supposed negative
		Negative	-	Supposed negative

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10. 引用文献（最終アクセス日：2019年11月1日）
- 1) OECD Guidelines for the Testing of Chemicals, Section 4, Available at: http://www.oecd-ilibrary.org/environment/oecd-guidelines-for-the-testing-of-chemicals-section-4-health-effects_20745788
 - 2) Corrigendum to Regulation (EC) No 1907/2006 of the European Parliament and of the Council of 18 December 2006 concerning the Registration, Evaluation, Authorisation and Restriction of Chemicals (REACH), establishing
 - 3) REGULATION (EC) No 1223/2009 OF THE EUROPEAN PARLIAMENT AND OF THE COUNCIL of 30 November 2009 on cosmetic products (recast) (Text with EEA relevance), Available at: https://ec.europa.eu/health/sites/health/files/endocrine_disruptors/docs/cosmetic_1223_2009_regulation_en.pdf
 - 4) Takahashi T, Kimura Y, Saito R, Nakajima Y, Ohmiya Y, Yamasaki K, Aiba S. (2011) An in vitro test to screen skin sensitizers using a stable THP-1-derived IL-8 reporter cell line, THP-G8. *Toxicol Sci.*; 124:359-69
 - 5) OECD (2017). Validation report for the international validation study on the IL-8 Luc assay as a test evaluating the skin sensitizing potential of chemicals conducted by the IL-8 Luc Assay. Series on Testing and Assessment No. 267, ENV/JM/MONO (2017)19. Organization for Economic Cooperation and Development, Paris. Available at: <https://www.oecd.org/env/ehs/testing/series-testing-assessment-publications-number.htm>
 - 6) Report of the peer review panel for the IL-8 Luciferase (IL-8 Luc) assay for in vitro skin sensitization. Series on testing and assessment. No.268 ENV/JM/MONO(2017)20. Organization for Economic Cooperation and Development, Paris. Available at: <http://www.oecd.org/env/ehs/testing/series-testing-assessment-publications-number.htm>
 - 7) OECD (2014) The Adverse Outcome Pathway for Skin Sensitisation Initiated by Covalent Binding to Proteins, OECD Publishing, Paris, Available at: <http://dx.doi.org/10.1787/9789264221444-en>
 - 8) Nukada Y, Miyazawa M, Kosaka N, Ito Y, Sakaguchi H, Nishiyama N. (2008) Production of IL-8 in THP-1 cells following contact allergen stimulation via mitogen-activated protein kinase activation or tumor necrosis factor-alpha production. *J Toxicol Sci.*, 33, 175–85.
 - 9) Weber FC1, Németh T, Csepregi JZ, Dudeck A, Roers A, Ozsvári B, Oswald E, Puskás LG, Jakob T, Mócsai A, Martin SF. (2015) Neutrophils are required for both the sensitization and elicitation phase of contact hypersensitivity, *J. Exp. Med.*, 212 (1), 15-22.
 - 10) Liu T, Zhang L, Joo D, Sun SC. (2017) NF- κ B signaling in inflammation, *Signal Transduct Target Ther.* 2, 17023.
 - 11) Kimura Y, Watanabe M, Suzuki N, Iwaki T, Yamakage K, Saito K, Nakajima Y, Fujimura C, Ohmiya Y, Omori T, Kojima H, Aiba S. (2018) The performance of an in vitro skin sensitisation test, IL-8 Luc assay (OECD442E), and the integrated approach with direct peptide reactive assay (DPRA), *J. Toxicol. Sci.*, 43 (12), 741-749.
 - 12) Ashikaga T, Sakaguchi H, Sono S, Kosaka N, Ishikawa M, Nukada Y, Miyazawa M, Ito Y, Nishiyama N, and Itagaki H. (2010) A comparative evaluation of in vitro skin sensitisation tests: the human cell-line activation test (h-CLAT) versus the local lymph node assay (LLNA). *Alternatives to laboratory animals : ATLA* 38:275-84.
 - 13) Guidance Document on the Reporting of Defined Approaches and Individual Information Sources to be Used within Integrated Approaches to Testing and Assessment (IATA) for Skin Sensitisation, Annex 1 and Annex 2. Available at: <http://www.oecd.org/chemicalsafety/risk-assessment/iata-integrated-approaches-to-testing-and-assessment.htm>

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Phase	Chemicals Repeat		No.	Chemical name	CAS No.	LLNA EC3	GHS Category	Lab A Judge	Lab B Judge	Lab C Judge
II b	5 coded	3 sets	1	p-Benzoquinone	106-51-4	0.099	1A	G(1/2)	G(1/2)	P(2/2)
								P(2/2)	P(2/2)	P(2/2)
								G(1/2)	P(2/2)	P(2/2)
			2	1-Thioglycerol	96-27-5	3.6	1B	P(4/4)	P(4/4)	P(4/4)
								P(2/3)	P(3/3)	P(3/3)
								P(3/3)	P(2/2)	P(3/3)
			3	Benzyl cinnamate	103-41-3	18.4	1B	P(4/4)	P(2/2)	P(4/4)
								P(3/3)	P(2/2)	P(3/3)
								P(4/4)	P(3/3)	P(3/3)
			4	Diethyl maleate	141-05-9	5.8	1B	N(0/3)	N(0/4)	N(0/3)
								N(0/3)	N(0/3)	N(0/4)
								N(0/3)	N(0/3)	N(1/4)
			5	2,4-Dichloronitrobenzene	611-06-3		NC	N(0/3)	N(0/3)	G(1/3)
								N(0/3)	G(0/2)	N(0/4)
								N(0/4)	G(1/3)	N(1/4)
II c	5 coded	3 sets	1	1,4-Phenylenediamine	106-50-3	0.0099	1A	P(2/2)	P(2/4)	N(1/4)
								P(2/4)	N(1/4)	N(1/4)
								P(3/4)	P(2/2)	P(2/3)
			2	2-Mercaptobenzothiazole	149-30-4	1.7	1A	P(2/2)	G(1/3)	P(2/3)
								P(2/2)	P(2/2)	P(2/2)
								P(2/2)	P(2/2)	P(2/2)
			3	Glycerol	56-81-5		NC	N(0/3)	N(0/4)	N(0/4)
								N(1/4)	N(0/3)	N(0/3)
								N(1/4)	N(1/4)	N(0/4)
			4	Isopropanol	67-63-0		NC	N(1/4)	P(2/2)	N(1/4)
								N(0/3)	P(2/2)	P(2/4)
								N(0/3)	P(2/4)	P(2/4)
			5	Methyl salicylate	119-36-8		NC	N(0/3)	N(0/3)	N(0/3)
								P(2/4)	N(0/4)	N(0/4)
								N(0/3)	N(1/4)	N(0/3)

III	20 coded	1 set	1	2,4-Dichloronitrobenzene	611-06-3		NC	N(1/4)	N(0/3)	N(0/3)
			2	DNCB(1-Chloro-2,4-dinitrobenzene)	97-00-7	0.06	1A	P(2/2)	P(2/2)	P(2/2)
			3	Glyoxal solution	107-22-2	0.8	1A	P(2/2)	G(1/3)	P(2/2)
			4	1,2-Benzisothiazol-3(2H)-one	2634-33-5	8	1B	P(2/2)	P(2/2)	P(2/2)
			5	Citral	5392-40-5	5	1B	P(2/2)	P(2/2)	P(2/2)
			6	4-Hydroxybenzoic acid	99-96-7		NC	N(0/4)	N(0/3)	N(0/4)
			7	Beryllium sulphate	7787-56-6	0.001	1A	N(1/4)	N(0/3)	N(0/3)
			8	Formaldehyde	50-00-0	0.61	1A	P(2/3)	P(2/2)	P(2/3)
			9	Chlorpromazine HCL	69-09-0	0.14	1A	P(2/2)	P(2/2)	P(2/2)
			10	Chloramine T	127-65-1	0.4	1A	P(2/3)	P(2/2)	P(2/2)
			11	Benzyl salicylate	118-58-1	2.9	1B	P(2/2)	P(2/2)	P(2/3)
			12	Dihydroeugenol	2785-87-7	6.8	1B	P(2/3)	P(2/2)	P(2/2)
			13	Nickel chloride	7718-54-9		1B	P(2/2)	P(2/2)	P(2/2)
			14	Imidazolidinyl urea	39236-46-9	24	1B	P(2/3)	P(2/3)	P(2/2)
			15	R(+)-Limonene	5989-27-5	69	1B	P(2/2)	N(0/4)	N(0/3)
			16	Methylmethacrylate	80-62-6	90	1B	N(0/3)	N(0/4)	N(0/4)
			17	Benzyl alcohol	100-51-6		NC	N(0/4)	N(0/3)	N(0/3)
			18	Dimethyl isophthalate	1459-93-4		NC	N(0/3)	N(0/4)	N(1/4)
			19	4-Aminobenzoic acid	150-13-0		NC	N(0/3)	N(0/3)	N(0/4)
			20	Xylene	1330-20-7	95.8	NC	N(0/4)	N(0/3)	N(0/3)
IV	5 coded	3 sets	1	Ethylene glycol dimethacrylate	97-90-5	28	1B	P(2/2)	P(2/4)	P(2/2)
								P(2/2)	P(2/2)	P(2/2)
								P(2/2)	P(2/2)	P(2/2)
			2	4-Nitrobenzylbromide	100-11-8	0.05	1A	P(2/2)	P(2/2)	P(2/2)
								P(2/2)	P(2/2)	P(2/2)
								P(2/2)	P(2/2)	P(2/2)
			3	4-Allylanisole	140-67-0	18	1B	P(2/2)	P(2/2)	P(2/2)
								P(2/2)	P(2/2)	P(2/2)
								P(2/2)	P(2/2)	P(2/2)
			4	Ethylendiamine	107-15-3	2.2	1B	P(2/4)	N(1/4)	N(0/3)
								P(2/3)	N(1/4)	N(1/4)
								P(2/2)	N(1/4)	N(1/4)
			5	Salicylic acid	69-72-7		NC	N(0/3)	N(1/4)	N(1/4)
								P(2/3)	N(1/4)	N(0/3)
								N(1/4)	N(0/3)	N(0/3)

P: positive, N: negative, G: inconclusive

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No.	CAS#	Chemical	GHS	LLNA Potency category	IL-8 Luc assay (DMSO)	IL-8 Luc assay (X-VIVO)	IL-8 Luc assay (Test Guideline protocol) (Op. 2)
1	93-91-4	1-Benzoylacetone (1-BA)	1A	Extreme	positive	positive	positive
2	1086-00-6	1-Chloromethylpyrene	1A	Extreme	positive	positive	positive
3	97-00-7	2,4-Dinitrochlorobenzene (DNCEB)	1A	Extreme	positive	positive	positive
4	100-11-8	4-Nitrobenzyl bromide (4-NBB)	1A	Extreme	positive	positive	positive
5	20048-27-5	Bandrowski's base	1A	Extreme	positive	positive	positive
6	7787-56-6	Beryllium sulfate tetrahydrate	1A	Extreme	negative	positive	positive
7	70-34-8	Dinitrofluorobenzene	1A	Extreme	positive	positive	positive
8	886-38-4	Diphenylcyclopropenone (DPCP)	1A	Extreme	positive	positive	positive
9	15646-46-5	Oxazolone	1A	Extreme	negative	positive	positive
10	106-51-4	p-Benzoquinone	1A	Extreme	positive	positive	positive
11	7778-50-9	Potassium dichromate	1A	Extreme	positive	positive	positive
12	1154-59-2	Tetrachlorosalicylanilide	1A	Extreme	positive	positive	positive
13	123-31-9	1,4-Dihydroquinone	1A	Strong	positive	positive	positive
14	106-50-3	1,4-Phenylenediamine(PPD)	1A	Strong	negative	positive	positive
15	95-55-6	2-Aminophenol	1A	Strong	positive	positive	positive
16	5307-14-2	2-Nitro-1,4-phenylenediamine	1A	Strong	positive	negative	inconclusive
17	94-36-0	Benzoyl peroxide	1A	Strong	positive	negative	inconclusive
18	100-39-0	Benzyl bromide	1A	Strong	positive	positive	positive
19	127-65-1	Chloramine-T	1A	Strong	positive	positive	positive
20	69-09-0	Chlorpromazine hydrochloride	1A	Strong	positive	positive	positive
21	7646-79-9	Cobalt chloride hexahydrate	1A	Strong	positive	positive	positive
22	50-00-0	Formaldehyde	1A	Strong	positive	positive	positive
23	111-30-8	Glutaraldehyde	1A	Strong	positive	positive	positive
24	6259-76-3	Hexyl salicylate	1A	Strong	positive	negative	inconclusive
25	55406-53-6	Iodopropynyl butylcarbamate	1A	Strong	positive	positive	positive
26	1166-52-5	Lauryl gallate	1A	Strong	positive	positive	positive
27	108-31-6	Maleic anhydride	1A	Strong	negative ^{out of AD}	negative ^{out of AD}	inconclusive, negative ^{out of AD}
28	35691-65-7	Methyldibromoglutaronitrile (MDGN)	1A	Strong	positive	positive	positive
29	85-44-9	Phthalic anhydride	1A	Strong	negative ^{out of AD}	negative ^{out of AD}	inconclusive, negative ^{out of AD}
30	121-79-9	Propyl gallate	1A	Strong	positive	positive	positive
31	90-15-3	1-Naphthol	1A	Moderate	positive	positive	positive
32	579-07-7	1-Phenyl-1,2-propanedione	1A	Moderate	positive	positive	positive
33	96-27-5	1-Thioglycerol(3-Mercapto-1,2-propanediol)	1B	Moderate	positive	positive	positive
34	2634-33-5	1,2-Benzisothiazolin-3-one	1B	Moderate	positive	positive	positive

No.	CAS#	Chemical	GHS	LLNA Potency category	IL-8 Luc assay (DMSO)	IL-8 Luc assay (X-VIVO)	IL-8 Luc assay (Test Guideline protocol) (Op. 2)
35	3344-77-2	12-Bromo-1-dodecanol	1B	Moderate	positive	positive	positive
36	818-61-1	2-Hydroxyethyl acrylate	1A	Moderate	negative	negative	negative
37	149-30-4	2-Mercaptobenzothiazole (2-MBT)	1A	Moderate	positive	positive	positive
38	93-51-6	2-Methoxy-4-methylphenol	1B	Moderate	positive	positive	positive
39	93-53-8	2-Phenylpropionaldehyde	1B	Moderate	positive	positive	positive
40	4313-03-5	2,4-Heptadienal	1B	Moderate	positive	positive	positive
41	591-27-5	3-Aminophenol	1B	Moderate	positive	positive	positive
42	109-55-7	3-Dimethylaminopropylamine	1B	Moderate	positive	positive	positive
43	17369-59-4	3-Propylidene-phthalide	1B	Moderate	positive	positive	positive
44	119-84-6	3,4-Dihydrocoumarin	1B	Moderate	negative	positive	positive
45	106-47-8	4-Chloroaniline	1B	Moderate	negative	positive	positive
46	21834-92-4	5-Methyl-2-phenyl-2-hexenal	1B	Moderate	positive	positive	positive
47	101-39-3	a-Methylcinnamic aldehyde	1B	Moderate	positive	positive	positive
48	118-58-1	Benzyl salicylate	1B	Moderate	positive	positive	positive
49	122-57-6	Benzylideneacetone	1B	Moderate	positive	positive	positive
50	104-55-2	Cinnamic aldehyde (Cinnamal)	1B	Moderate	positive	positive	positive
51	23593-75-1	Clotrimazole	1B	Moderate	positive	negative	inconclusive
52	141-05-9	Diethyl maleate	1B	Moderate	negative	positive	positive
53	64-67-5	Diethyl sulfate	1B	Moderate	negative	positive	positive
54	111-40-0	Diethylenetriamine	1B	Moderate	positive	positive	positive
55	2785-87-7	Dihydroeugenol (2-Methoxy-4-propylphenol)	1B	Moderate	positive	positive	Positive
56	107-15-3	Ethylenediamine	1B	Moderate	negative	positive	Positive
57	107-22-2	Glyoxal	1A	Moderate	negative	positive	Positive
58	97-54-1	Isoeugenol	1A	Moderate	positive	positive	Positive
59	111-80-8	Methyl-2-nonynoate	1B	Moderate	negative	positive	Positive
60	2682-20-4	Methylisothiazolinone, 2-Methyl-2H-Isothiazol-3-one	1A	Moderate	positive	positive	Positive
61	10101-97-0	Nickel sulphate hexahydrate	1B	Moderate	positive	positive	Positive
62	7718-54-9	Nickel (II) chloride		Moderate	positive	positive	Positive
63	18031-40-8	Perillaldehyde	1B	Moderate	positive	positive	Positive
64	122-78-1	Phenylacetaldehyde	1B	Moderate	positive	positive	Positive
65	108-46-3	Resorcinol	1B	Moderate	positive	positive	Positive
66	133-37-9	Tartaric acid	1B	Moderate	positive	positive	Positive
67	137-26-8	Tetramethylthiuramdisulphide (TMTD)	1B	Moderate	positive	positive	Positive
68	6728-26-3	trans-2-Hexenal	1B	Moderate	positive	positive	Positive
69	552-30-7	Trimellitic anhydride	1B	Moderate	negative	positive	Positive

No.	CAS#	Chemical	GHS	LLNA Potency category	IL-8 Luc assay (DMSO)	IL-8 Luc assay (X-VIVO)	IL-8 Luc assay (Test Guideline protocol) (Op. 2)
70	112-45-8	Undec-10-enal	1B	Moderate	positive	positive	Positive
71	5989-27-5	(R)-(+)-Limonene	1B	Weak	negative	positive	Positive
72	111-25-1	1-Bromohexane	1B	Weak	negative	negative	inconclusive
73	103-11-7	2-Ethylhexyl acrylate	1B	Weak	positive	positive	Positive
74	431-03-8	2,3-Butanedione	1B	Weak	positive	positive	Positive
75	140-67-0	4-Allylanisole	1B	Weak	negative	positive	Positive
76	13706-86-0	5-Methyl-2,3-hexanedione	1B	Weak	negative	positive	Positive
77	514-10-3	Abietic acid	1B	Weak	positive	positive	Positive
78	122-40-7	Amyl cinnamic aldehyde	1B	Weak	positive	positive	Positive
79	62-53-3	Aniline	1B	Weak	negative	negative	Negative
80	94-09-7	Benzocaine	1B	Weak	negative	negative	inconclusive
81	103-41-3	Benzyl cinnamate	1B	Weak	positive	positive	Positive
82	1565-94-2	Bis-GMA	1B	Weak	positive	positive	Positive
83	2426-08-6	Butyl glycidyl ether	1B	Weak	positive	positive	Positive
84	104-54-1	Cinnamic alcohol	1B	Weak	positive	positive	Positive
85	5392-40-5	Citral	1B	Weak	positive	positive	Positive
86	103-95-7	Cyclamen aldehyde	1B	Weak	positive	positive	Positive
87	97-90-5	Ethylene glycol dimethacrylate	1B	Weak	negative	positive	Positive
88	97-53-0	Eugenol	1B	Weak	positive	positive	Positive
89	19317-11-4	Farnesal	1B	Weak	positive	positive	Positive
90	106-24-1	Geraniol	1B	Weak	positive	positive	Positive
91	101-86-0	Hexyl cinnamic aldehyde	1B	Weak	positive	positive	Positive
92	107-75-5	Hydroxycitronellal	1B	Weak	positive	positive	Positive
93	39236-46-9	Imidazolidinyl urea	1B	Weak	positive	positive	Positive
94	80-54-6	Lilial	1B	Weak	positive	positive	Positive
95	78-70-6	Linalool	1B	Weak	positive	positive	Positive
96	31906-04-4	Lylal	1B	Weak	positive	positive	Positive
97	80-62-6	Methyl methacrylate	1B	Weak	negative	negative	inconclusive
98	613-29-6	N,N-Dibutylaniline	1B	Weak	positive	positive	Positive
99	112-05-0	Nonanoic acid	1B	Weak	positive	positive	Positive
100	144-62-7	Oxalic acid	1B	Weak	positive	positive	Positive
101	61-33-6	Penicillin G	1B	Weak	negative	negative	Negative
102	87-86-5	Pentachlorophenol	1B	Weak	positive	positive	Positive
103	2050-14-8	Phenol, 2,2-azobis-	1B	Weak	positive	positive	Positive
104	93-99-2	Phenyl benzoate	1B	Weak	positive	negative	inconclusive
105	110-86-1	Pyridine	1B	Weak	negative	negative	Negative

No.	CAS#	Chemical	GHS	LLNA Potency category	IL-8 Luc assay (DMSO)	IL-8 Luc assay (X-VIVO)	IL-8 Luc assay (Test Guideline protocol) (Op. 2)
106	1330-20-7	Xylene	1B	Weak	negative	negative	inconclusive
107	127-51-5	α -iso-Methylionone	1B	Weak	positive	positive	Positive
108	109-65-9	1-Bromobutane	Non	negative	negative	negative	inconclusive
109	71-36-3	1-Butanol	Non	negative	negative	negative	Negative
110	874-23-7	2-Acetylcyclohexanone	Non	negative	positive	positive	Positive
111	923-26-2	2-Hydroxypropyl methacrylate	Non	negative	negative	negative	Negative
112	611-06-3	2,4-Dichloro-1-nitrobenzene	Non	negative	negative	negative	inconclusive
113	99-96-7	4-Hydroxybenzoic acid	Non	negative	negative	negative	inconclusive
114	92-48-8	6-Methylcoumarin	Non	negative	negative	negative	inconclusive
115	100-06-1	Acetoanisole (4-Methoxyacetophenone)	Non	negative	negative	negative	inconclusive
116	8001-54-5	Benzalkonium chloride	Non	negative	positive ^{out} of AD	positive ^{out} of AD	positive ^{out} of AD
117	65-85-0	Benzoic acid	Non	negative	negative	negative	inconclusive
118	100-51-6	Benzyl alcohol	Non	negative	negative	positive	Positive
119	108-90-7	Chlorobenzene	Non	negative	negative	negative	inconclusive
120	9004-54-0	Dextran	Non	negative	negative	negative	Negative
121	84-66-2	Diethyl phthalate	Non	negative	positive	positive	Positive
122	68-12-2	Dimethyl formamide	Non	negative	negative	negative	Negative
123	1459-93-4	Dimethyl isophthalate	Non	negative	negative	negative	inconclusive
124	94-02-0	Ethyl benzoylacetate	Non	negative	positive	negative	inconclusive
125	121-32-4	Ethyl vanilline	Non	negative	positive	positive	positive
126	492-94-4	Furil	Non	negative	positive	negative	inconclusive
127	56-81-5	Glycerol	Non	negative	negative	negative	negative
128	57-09-0	Hexadecyltrimethylammonium bromide	Non	negative	positive ^{out} of AD	positive ^{out} of AD	positive ^{out} of AD
129	67-63-0	Isopropanol	Non	negative	negative	negative	negative
130	50-21-5	Lactic acid	Non	negative	negative	negative	negative
131	119-36-8	Methyl salicylate	Non	negative	negative	negative	inconclusive
132	124-07-2	Octanoic acid (Caprylic acid)	Non	negative	positive ^{out} of AD	positive ^{out} of AD	positive ^{out} of AD
133	150-13-0	p-Aminobenzoic acid	Non	negative	negative	negative	inconclusive
134	94-13-3	Propyl paraben	Non	negative	negative	positive	positive
135	57-55-6	Propylene glycol	Non	negative	negative	negative	negative
136	81-07-2	Saccharin	Non	negative	negative	positive	positive
137	69-72-7	Salicylic acid	Non	negative	negative	negative	inconclusive
138	151-21-3	Sodium lauryl sulfate	Non	negative	positive ^{out} of AD	positive ^{out} of AD	positive ^{out} of AD
139	3810-74-0	Streptomycin sulfate	Non	negative	negative	negative	negative
140	63-74-1	Sulfanilamide	Non	negative	negative	negative	inconclusive

No.	CAS#	Chemical	GHS	LLNA Potency category	IL-8 Luc assay (DMSO)	IL-8 Luc assay (X-VIVO)	IL-8 Luc assay (Test Guideline protocol) (Op. 2)
141	9005-65-6	Tween-80	Non	negative	positive ^{out} of AD	positive ^{out} of AD	positive ^{out} of AD
142	121-33-5	Vanillin	Non	negative	positive	positive	positive
143	7733-02-0	Zinc sulfate	Non	negative	positive	positive	positive

out of AD: out of applicability domain, inconclusive: Should not be considered

*KEY EVENT-BASED TEST GUIDELINE***IN VITRO SKIN SENSITISATION ASSAYS ADDRESSING THE KEY
EVENT ON ACTIVATION OF DENDRITIC CELLS ON THE
ADVERSE OUTCOME PATHWAY FOR SKIN SENSITISATION****GENERAL INTRODUCTION***Activation of dendritic cells Key Event based Test Guideline*

1. A skin sensitiser refers to a substance that will lead to an allergic response following skin contact as defined by the United Nations Globally Harmonized System of Classification and Labelling of Chemicals (UN GHS) (1). There is general agreement on the key biological events underlying skin sensitisation. The current knowledge of the chemical and biological mechanisms associated with skin sensitisation has been summarised as an Adverse Outcome Pathway (AOP) (2), starting with the molecular initiating event through intermediate events to the adverse effect, namely allergic contact dermatitis. In this instance, the molecular initiating event (i.e. the first key event) is the covalent binding of electrophilic substances to nucleophilic centres in skin proteins. The second key event in this AOP takes place in the keratinocytes and includes inflammatory responses as well as changes in gene expression associated with specific cell signalling pathways such as the antioxidant/electrophile response element (ARE)-dependent pathways. The third key event is the activation of dendritic cells (DC), typically assessed by expression of specific cell surface markers, chemokines and cytokines. The fourth key event is T-cell activation and proliferation, which is indirectly assessed in the murine Local Lymph Node Assay (LLNA) (3).

2. This Test Guideline (TG) describes *in vitro* assays that address mechanisms described under the Key Event on activation of dendritic cells of the AOP for skin sensitisation (2). The TG comprises test methods to be used for supporting the discrimination between skin sensitisers and non-sensitisers in accordance with the UN GHS (1).

The test methods described in this TG are:

- Human Cell Line Activation test (h-CLAT)
- U937 cell line activation Test (U-SENSTM)
- Interleukin-8 Reporter Gene Assay (IL-8 Luc assay)

© OECD, (2018)

You are free to use this material subject to the terms and conditions available at <http://www.oecd.org/termsandconditions/>.

In accordance with the decision of the Council on a delegation of authority to amend Annex I of the decision of the council on the Mutual Acceptance of Data in the assessment of chemicals [C(2018)49], this Guideline was amended by the OECD's Joint Meeting of the Chemicals Committee and the Working Party on Chemicals, Pesticides and Biotechnology by written procedure on 25 June 2018.

3. The test methods included in this Test Guideline may differ in relation to the procedure used to generate the data and the readouts measured but can be used indiscriminately to address countries' requirements for test results on the Key Event on activation of dendritic cells of the AOP for skin sensitisation while benefiting from the Mutual Acceptance of Data.

Background and principles of the test methods included in the Key Event based Test Guideline

4. The assessment of skin sensitisation has typically involved the use of laboratory animals. The classical methods that use guinea-pigs, the Guinea Pig Maximisation Test (GPMT) of Magnusson and Kligman, and the Buehler Test (TG 406) (4), assess both the induction and elicitation phases of skin sensitisation. The murine tests, the LLNA (TG 429) (3) and its two non-radioactive modifications, LLNA: DA (TG 442 A) (5) and LLNA: BrdU-ELISA (TG 442 B) (6), all assess the induction response exclusively, and have also gained acceptance, since they provide an advantage over the guinea pig tests in terms of animal welfare together with an objective measurement of the induction phase of skin sensitisation.

5. Recently mechanistically-based *in chemico* and *in vitro* test methods addressing the first key event (OECD TG 442C; Direct Peptide Reactivity Assay (7)), and second key event (OECD TG 442D; ARE-Nrf2 Luciferase Test Method (8)) of the skin sensitisation AOP have been adopted for contributing to the evaluation of the skin sensitisation hazard potential of chemicals.

6. Test methods described in this TG either quantify the change in the expression of cell surface marker(s) associated with the process of activation of monocytes and DC following exposure to sensitisers (e.g. CD54, CD86) or the changes in IL-8 expression, a cytokine associated with the activation of DC. Skin sensitisers have been reported to induce the expression of cell membrane markers such as CD40, CD54, CD80, CD83, and CD86 in addition to induction of proinflammatory cytokines, such as IL-1 β and TNF- α , and several chemokines including IL-8 (CXCL8) and CCL3 (9) (10) (11) (12), associated with DC activation (2).

7. However, as DC activation represents only one key event of the skin sensitisation AOP (2) (13), information generated with test methods measuring markers of DC activation alone may not be sufficient to conclude on the presence or absence of skin sensitisation potential of chemicals. Therefore data generated with the test methods described in this Test Guideline are proposed to support the discrimination between skin sensitisers (i.e. UN GHS Category 1) and non-sensitisers when used within Integrated Approaches to Testing and Assessment (IATA), together with other relevant complementary information, e.g. derived from *in vitro* assays addressing other key events of the skin sensitisation AOP as well as non-testing methods, including read-across from chemical analogues (13). Examples of the use of data generated with these methods within Defined Approaches, i.e. approaches standardised both in relation to the set of information sources used and in the procedure applied to the data to derive predictions, have been published (13) and can be employed as useful elements within IATA.

8. The test methods described in this Test Guideline cannot be used on their own, neither to sub-categorise skin sensitisers into subcategories 1A and 1B as defined by UN GHS (1), for authorities implementing these two optional subcategories, nor to predict potency for safety assessment decisions. However, depending on the regulatory framework, positive results generated with these methods may be used on their own to classify a chemical into UN GHS category 1.

9. The term "test chemical" is used in this Test Guideline to refer to what is being tested¹ and is not related to the applicability of the test methods to the testing of mono-constituent substances, multi-constituent substances and/or mixtures. Limited information is currently available on the applicability of the test methods to multi-constituent substances/mixtures (14) (15). The test methods are nevertheless technically applicable to the testing of multi-constituent substances and mixtures. When considering testing of mixtures, difficult-to-test chemicals (e.g. unstable), or test chemicals not clearly within the applicability domain described in this Guideline, upfront consideration should be given to whether the results of such testing will yield results that are meaningful scientifically. Moreover, when testing multi-constituent substances or mixtures, consideration should be given to possible interference of cytotoxic constituents with the observed responses.

LITERATURE

1. United Nations UN (2015), Globally Harmonized System of Classification and Labelling of Chemicals (GHS). Sixth revised edition. New York & Geneva: United Nations Publications. ISBN: 978-92-1-117006-1.
Available at: [https://www.unece.org/trans/danger/publi/ghs/ghs_rev06/06files_e.html].
2. OECD (2012), The Adverse Outcome Pathway for Skin Sensitisation Initiated by Covalent Binding to Proteins. Part 1: Scientific Evidence. Series on Testing and Assessment No. 168. Available at: [http://www.oecd.org/officialdocuments/publicdisplaydocumentpdf/?cote=ENV/JM/MONO\(2012\)10/PART1&docLanguage=En](http://www.oecd.org/officialdocuments/publicdisplaydocumentpdf/?cote=ENV/JM/MONO(2012)10/PART1&docLanguage=En)
3. OECD (2010), The Local Lymph Node Assay. OECD Guideline for the Testing of Chemicals No. 429, Organisation for Economic Cooperation and Development, Paris. Available at: <http://www.oecd.org/env/testguidelines>
4. OECD (1992), Skin Sensitisation. OECD Guideline for the Testing of Chemicals No. 406, Organisation for Economic Cooperation and Development, Paris. Available at: <http://www.oecd.org/env/testguidelines>
5. OECD (2010), Skin Sensitisation: Local Lymph Node Assay: DA. OECD Guideline for the Testing of Chemicals No. 442A, Organisation for Economic Cooperation and Development, Paris. Available at: <http://www.oecd.org/env/testguidelines>
6. OECD (2010), Skin sensitisation: Local Lymph Node Assay: BrdU-ELISA. OECD Guideline for the Testing of Chemicals No. 442B, Organisation for Economic Cooperation and Development, Paris. Available at: <http://www.oecd.org/env/testguidelines>
7. OECD (2015), OECD Guideline for the Testing of Chemicals No. 442C: In Chemico Skin Sensitisation: Direct Peptide Reactivity Assay (DPRA). Paris, France: Organisation for Economic Cooperation and Development. Available at: <http://www.oecd.org/env/testguidelines>

¹ In June 2013, the Joint Meeting agreed that where possible, a more consistent use of the term "test chemical" describing what is being tested should be applied in new and updated Test Guidelines.

8. OECD (2015), OECD Guideline for the Testing of Chemicals No. 442D: In Vitro Skin Sensitisation: ARE-Nrf2 Luciferase Test Method. Paris, France: Organisation for Economic Cooperation and Development. Available at: <http://www.oecd.org/env/testguidelines>
9. Steinman RM. 1991. The dendritic cell system and its role in immunogenicity. *Annu Rev Immunol* 9:271-96.
10. Caux C, Vanbervliet B, Massacrier C, Azuma M, Okumura K, Lanier LL, and Banchereau J. 1994. B70/B7-2 is identical to CD86 and is the major functional ligand for CD28 expressed on human dendritic cells. *J Exp Med* 180:1841-7.
11. Aiba S, Terunuma A, Manome H, and Tagami H. 1997. Dendritic cells differently respond to haptens and irritants by their production of cytokines and expression of co-stimulatory molecules. *Eur J Immunol* 27:3031-8.
12. Aiba S, Manome H, Nakagawa S, Mollah ZU, Mizuashi M, Ohtani T, Yoshino Y, and Tagami H. 2003. p38 mitogen-activated protein kinase and extracellular signal-regulated kinases play distinct roles in the activation of dendritic cells by two representative haptens, NiCl₂ and DNCB. *J Invest Dermatol* 120:390-8.
13. OECD (2016). Series on Testing & Assessment No. 256: Guidance Document On The Reporting Of Defined Approaches And Individual Information Sources To Be Used Within Integrated Approaches To Testing And Assessment (IATA) For Skin Sensitisation, Annex 1 and Annex 2. [ENV/JM/HA\(2016\)29](#). Organisation for Economic Cooperation and Development, Paris. Available at: [<https://community.oecd.org/community/iatass>].
14. Ashikaga T, Sakaguchi H, Sono S, Kosaka N, Ishikawa M, Nukada Y, Miyazawa M, Ito Y, Nishiyama N, Itagaki H. (2010), A comparative evaluation of in vitro skin sensitisation tests: the human cell-line activation test (h-CLAT) versus the local lymph node assay (LLNA). *Altern. Lab. Anim.* 38, 275-284.
15. Piroird, C., Ovigne, J.M., Rousset, F., Martinozzi-Teissier, S., Gomes, C., Cotovio, J., Alépée, N. (2015). The Myeloid U937 Skin Sensitization Test (U-SENS) addresses the activation of dendritic cell event in the adverse outcome pathway for skin sensitization. *Toxicol. In Vitro* 29, 901-916.

APPENDIX II

PROFICIENCY SUBSTANCES

Prior to routine use of the test method described in this Annex to Test Guideline 442E, laboratories should demonstrate technical proficiency by correctly obtaining the expected U-SENS™ prediction for the 10 substances recommended in Table 1 and by obtaining CV70 and EC150 values that fall within the respective reference range for at least 8 out of the 10 proficiency substances. Proficiency substances were selected to represent the range of responses for skin sensitisation hazards. Other selection criteria were that the substances are commercially available, and that high-quality *in vivo* reference data as well as high quality *in vitro* data generated with the U-SENS™ method are available. Also, published reference data are available for the U-SENS™ method (1) (8).

Table 1: Recommended substances for demonstrating technical proficiency with the U-SENS™ method

Proficiency substances	CASRN	Physical state	<i>In vivo</i> prediction ¹	U-SENS™ Solvent/ Vehicle	U-SENS™ CV70 Reference Range in µg/mL ²	U-SENS™ EC150 Reference Range in µg/mL ²
4-Phenylenediamine	106-50-3	Solid	Sensitiser (strong)	Complete medium ³	<30	Positive (≤10)
Picryl sulfonic acid	2508-19-2	Liquid	Sensitizer (strong)	Complete medium	>50	Positive (≤50)
Diethyl maleate	141-05-9	Liquid	Sensitiser (moderate)	DMSO	10-100	Positive (≤20)
Resorcinol	108-46-3	Solid	Sensitiser (moderate)	Complete medium	>100	Positive (≤50)
Cinnamic alcohol	104-54-1	Solid	Sensitiser (weak)	DMSO	>100	Positive (10-100)
4-Allylanisole	140-67-0	Liquid	Sensitiser (weak)	DMSO	>100	Positive (≤200)
Saccharin	81-07-2	Solid	Non-sensitiser	DMSO	>200	Negative (>200)
Glycerol	56-81-5	Liquid	Non-sensitiser	Complete medium	>200	Negative (>200)
Lactic acid	50-21-5	Liquid	Non-sensitiser	Complete medium	>200	Negative (>200)
Salicylic acid	69-72-7	Solid	Non-sensitiser	DMSO	>200	Negative (>200)

Abbreviations: CAS RN = Chemical Abstracts Service Registry Number

¹ The *in vivo* hazard and (potency) prediction is based on LLNA data (1) (8). The *in vivo* potency is derived using the criteria proposed by ECETOC (17).

² Based on historical observed values (1) (8).

³ Complete medium: RPMI-1640 medium supplemented with 10% foetal calf serum, 2 mM L-glutamine, 100 units/mL penicillin and 100 µg/mL streptomycin (8).

ANNEX III: IN VITRO SKIN SENSITISATION: IL-8 LUC ASSAY

INITIAL CONSIDERATIONS AND LIMITATIONS

1. In contrast to assays analysing the expression of cell surface markers, the IL-8-Luc assay quantifies changes in IL-8 expression, a cytokine associated with the activation of dendritic cells (DC). In the THP-1-derived IL-8 reporter cell line (THP-G8, established from the human acute monocytic leukemia cell line THP-1), IL-8 expression is measured following exposure to sensitisers (1). The expression of luciferase is then used to aid discrimination between skin sensitisers and non-sensitisers.

2. The IL-8 Luc method has been evaluated in a validation study (2) conducted by the Japanese Centre for the Validation of Alternatives Methods (JaCVAM), the **Ministry of Economy, Trade and Industry** (METI), and the Japanese Society **for Alternatives to Animal Experiments** (JSAAE) and subsequently subjected to independent peer review (3) under the auspices of JaCVAM and the Ministry of Health, Labour and Welfare (MHLW) with the support of the **International Cooperation on Alternative Test Methods** (ICATM). Considering all available evidence and input from regulators and stakeholders, the IL-8 Luc assay is considered useful as part of IATA to discriminate sensitisers from non-sensitisers for the purpose of hazard classification and labelling. Examples of the use of IL-8 Luc assay data in combination with other information are reported in the literature (4) (5) (6).

3. The IL-8 Luc assay proved to be transferable to laboratories experienced in cell culture and luciferase measurement. Within and between laboratory reproducibilities were 87.7% and 87.5%, respectively (2). Data generated in the validation study (2) and other published work (1) (6) show that versus the LLNA, the IL-8 Luc assay judged 118 out of 143 chemicals as positive or negative and judged 25 chemicals as inconclusive and the accuracy of the IL-8 Luc assay in distinguishing skin sensitisers (UN GHS Cat. 1) from non-sensitisers (UN GHS No Cat.) is 86% (101/118) with a sensitivity of 96% (92/96) and specificity of 41% (9/22). Excluding substances outside the applicability domain described below (paragraph 5), the IL-8 Luc assay judged 113 out of 136 chemicals as positive or negative and judged 23 chemicals as inconclusive and the accuracy of the IL-8 Luc assay is 89% (101/113) with sensitivity of 96% (92/96) and specificity of 53% (9/17). Using human data cited in Urbisch et al. (7), the IL-8 Luc assay judged 76 out of 90 chemicals as positive or negative and judged 14 chemicals as inconclusive and the accuracy is 80% (61/76), sensitivity is 93% (54/58) and specificity is 39% (7/18). Excluding substances outside the applicability domain, the IL-8 Luc assay judged 71 out of 84 chemicals as positive or negative and judged 13 chemicals as inconclusive and the accuracy is 86% (61/71) with sensitivity of 93% (54/58) and specificity of 54% (7/13). False negative predictions with the IL-8 Luc assay are more likely to occur with chemicals showing low/moderate skin sensitisation potency (UN GHS subcategory 1B) than those with high potency (UN GHS subcategory 1A) (6). Together, the information supports a role for the IL-8 Luc assay in the identification of skin sensitisation hazards. The accuracy given for the IL-8 Luc assay as a standalone test method is only for guidance, as the method should be considered in combination with other sources of information in the context of an IATA and in accordance with the provisions of paragraphs 7 and 8 in the General introduction. Furthermore, when evaluating non-animal methods for skin sensitisation, it should be remembered that the LLNA and other animal tests may not fully reflect the situation in humans.

4. On the basis of the data currently available, the IL-8 Luc assay was shown to be applicable to test chemicals covering a variety of organic functional groups, reaction mechanisms, skin sensitisation potency (as determined in *in vivo* studies) and physicochemical properties (2) (6).

5. Although the IL-8 Luc assay uses X-VIVO™ 15 as a solvent, it correctly evaluated chemicals with a Log K_{ow} >3.5 and those with a water solubility of around 100 µg/ mL as calculated by EPI Suite™ and its performance to detect sensitisers with poor water solubility is better than that of the IL-8 Luc assay using dimethyl sulfoxide (DMSO) as a solvent (2). However, negative results for test chemicals that are not dissolved at 20 mg/ml may produce false negative results due to their inability to dissolve in X-VIVO™ 15. Therefore, negative results for these chemicals should not be considered. A high false negative rate for anhydrides was seen in the validation study. Furthermore, because of the limited metabolic capability of the cell line (8) and the experimental conditions, pro-haptens (substances requiring metabolic activation) and pre-haptens (substances activated by air oxidation) might give negative results in the assay. However, although negative results for suspected pre/prohaptens should be interpreted with caution, the IL-8 Luc assay correctly judged 11 out of 11 pre-haptens, 6/6 pro-haptens, and 6/8 pre/prohaptens in the IL-8 Luc assay data set (2). Based on the recent comprehensive review on three non-animal methods (the DPRA, the KeratinoSens™ and the h-CLAT) to detect pre and prohaptens (9), and based on the fact that THP-G8 cells used in the IL-8 Luc assay is a cell line derived from THP-1 that is used in the h-CLAT, the IL-8 Luc assay may also contribute to increase the sensitivity of non-animal methods to detect pre and pro-haptens in the combination of other methods. Surfactants tested so far gave (false) positive results irrespective of their type (e.g. cationic, anionic or on-ionic). Finally, chemicals that interfere with luciferase can confound its activity/measurement, causing apparent inhibition or increased luminescence (10). For example, phytoestrogen concentrations higher than 1 µM were reported to interfere with luminescence signals in other luciferase-based reporter gene assays due to over-activation of the luciferase reporter gene. Consequently, luciferase expression obtained at high concentrations of phytoestrogens or compounds suspected of producing phytoestrogen-like activation of the luciferase reporter gene needs to be examined carefully (11). Based on the above, surfactants, anhydrides and chemicals interfering with luciferase are outside the applicability domain of this assay. In cases where there is evidence demonstrating the non-applicability of the IL-8 Luc assay to other specific categories of test chemicals, the method should not be used for those specific categories.

6. As described above, the IL-8 Luc assay supports discrimination of skin sensitisers from non-sensitisers. Further work, preferably based on human data, is required to determine whether IL-8 Luc results can contribute to potency assessment when considered in combination with other information sources.

7. Definitions are provided in Appendix I.

PRINCIPLE OF THE TEST

8. The IL-8 Luc assay makes use of a human monocytic leukemia cell line THP-1 that was obtained from the American Type Culture Collection (Manassas, VA, USA). Using this cell line, the Dept. of Dermatology, Tohoku University School of Medicine, established a THP-1-derived IL-8 reporter cell line, THP-G8, that harbours the Stable Luciferase Orange (SLO) and Stable Luciferase Red (SLR) luciferase genes under the control of the IL-8 and **glyceraldehyde 3-phosphate dehydrogenase (GAPDH)** promoters, respectively (1). This allows quantitative measurement of luciferase gene induction by detecting luminescence from well-established light producing luciferase substrates as an indicator of the activity of the IL-8 and GAPDH in cells following exposure to sensitising chemicals.

9. The dual-colour assay system comprises an orange-emitting luciferase (SLO; $\lambda_{\text{max}} = 580 \text{ nm}$) (12) for the gene expression of the IL-8 promoter as well as a red-emitting luciferase (SLR; $\lambda_{\text{max}} = 630 \text{ nm}$) (13) for the gene expression of the internal control promoter, GAPDH. The two luciferases emit different colours upon reacting with firefly D-luciferin and their luminescence is measured simultaneously in a one-step reaction by dividing the emission from the assay mixture using an optical filter (14) (Appendix II).

10. THP-G8 cells are treated for 16 hours with the test chemical, after which SLO luciferase activity (SLO-LA) reflecting IL-8 promoter activity and SLR luciferase activity (SLR-LA) reflecting GAPDH promoter activity are measured. To make the abbreviations easy to understand, SLO-LA and SLR-LA are designated as IL8LA and GAPLA, respectively. Table 1 gives a description of the terms associated with luciferase activity in the IL-8 Luc assay. The measured values are used to calculate the normalised IL8LA (nIL8LA), which is the ratio of IL8LA to GAPLA; the induction of nIL8LA (Ind-IL8LA), which is the ratio of the arithmetic means of quadruple-measured values of the nIL8LA of THP-G8 cells treated with a test chemical and the values of the nIL8LA of untreated THP-G8 cells; and the inhibition of GAPLA (Inh-GAPLA), which is the ratio of the arithmetic means of quadruple-measured values of the GAPLA of THP-G8 cells treated with a test chemical and the values of the GAPLA of untreated THP-G8 cells, and used as an indicator for cytotoxicity.

Table 1. Description of terms associated with the luciferase activity in the IL-8 Luc assay

Abbreviations	Definition
GAPLA	SLR luciferase activity reflecting GAPDH promoter activity
IL8LA	SLO luciferase activity reflecting IL-8 promoter activity
nIL8LA	IL8LA / GAPLA
Ind-IL8LA	nIL8LA of THP-G8 cells treated with chemicals / nIL8LA of untreated cells
Inh-GAPLA	GAPLA of THP-G8 treated with chemicals / GAPLA of untreated cells
CV05	The lowest concentration of the chemical at which Inh-GAPLA becomes < 0.05.

11. Performance standards (PS) (15) are available to facilitate the validation of modified *in vitro* IL-8 luciferase test methods similar to the IL-8 Luc assay and allow for timely amendment of this Test Guideline for their inclusion. Mutual Acceptance of Data (MAD) will only be guaranteed for test methods validated according to the PS, if these test methods have been reviewed and included in this Test Guideline by the OECD (16).

DEMONSTRATION OF PROFICIENCY

12. Prior to routine use of the test method described in this Annex to Test Guideline 442E, laboratories should demonstrate technical proficiency, using the 9 Proficiency Substances listed in Appendix III in compliance with the Good *in vitro* Method Practices (17). Moreover, test method users should maintain a historical database of data generated with the reactivity checks (see paragraph 15) and with the positive and solvent/vehicle controls (see paragraphs 21-24), and use these data to confirm the reproducibility of the test method in their laboratory is maintained over time.

PROCEDURE

13. The Standard Operating Procedure (SOP) for the IL-8 Luc assay is available and should be employed when performing the test (18). Laboratories willing to perform the test can obtain the recombinant THP-G8 cell line from GPC Lab. Co. Ltd., Tottori, Japan, upon signing a Material Transfer Agreement (MTA) in line with the conditions of the OECD template. The following paragraphs provide a description of the main components and procedures of the assay.

Preparation of cells

14. The THP-G8 cell line from GPC Lab. Co. Ltd., Tottori, Japan, should be used for performing the IL-8 Luc assay (see paragraphs 8 and 13). On receipt, cells are propagated (2-4 passages) and stored frozen as a homogeneous stock. Cells from this stock can be propagated up to a maximum of 12 passages or a maximum of 6 weeks. The medium used for propagation is the RPMI-1640 culture medium containing 10% foetal bovine serum (FBS), antibiotic/antimycotic solution (100U/mL of penicillin G, 100µg/mL of streptomycin and 0.25µg/mL of amphotericin B in 0.85% saline) (e.g. GIBCO Cat#15240-062), 0.15µg/mL Puromycin (e.g. CAS:58-58-2) and 300µg/mL G418 (e.g. CAS:108321-42-2).

15. Prior to use for testing, the cells should be qualified by conducting a reactivity check. This check should be performed 1-2 weeks or 2-4 passages after thawing, using the positive control, 4-nitrobenzyl bromide (4-NBB) (CAS:100-11-8, ≥ 99% purity) and the negative control, lactic acid (LA) (CAS:50-21-5, ≥85% purity). 4-NBB should produce a positive response to Ind-IL8LA (≥1.4), while LA should produce a negative response to Ind-IL8LA (<1.4). Only cells that pass the reactivity check are used for the assay. The check should be performed according to the procedures described in paragraphs 22-24.

16. For testing, THP-G8 cells are seeded at a density of 2 to 5 × 10⁵ cells/mL, and pre-cultured in culture flasks for 48 to 96 hours. On the day of the test, cells harvested from the culture flask are washed with RPMI-1640 containing 10% FBS without any antibiotics, and then, resuspended with RPMI-1640 containing 10% FBS without any antibiotics at 1 × 10⁶ cells/mL. Then, cells are distributed into a 96-well flat-bottom black plate (e.g. Costar Cat#3603) with 50µL (5 × 10⁴ cells/well).

Preparation of the test chemical and control substances

17. The test chemical and control substances are prepared on the day of testing. For the IL-8 Luc assay, test chemicals are dissolved in X-VIVO™ 15, a commercially available serum-free medium (Lonza, 04-418Q), to the final concentration of 20 mg/mL. X-VIVO™ 15 is added to 20 mg of test chemical (regardless of the chemical's solubility) in a microcentrifuge tube and brought to a volume of 1mL and then vortexed vigorously and shaken on a rotor at a maximum speed of 8 rpm for 30 min at an ambient temperature of about 20°C. Furthermore, if solid chemicals are still insoluble, the tube is sonicated until the chemical is dissolved completely or stably dispersed. For test chemicals soluble in X-VIVO™ 15, the solution is diluted by a factor of 5 with X-VIVO™ 15 and used as an X-VIVO™ 15 stock solution of the test chemical (4 mg/mL). For test chemicals not soluble in X-VIVO™ 15, the mixture is rotated again for at least 30 min, then centrifuged at 15,000 rpm (≈20,000g) for 5 min; the resulting supernatant is used as an X-VIVO™ 15 stock solution of the test chemical. A scientific rationale should be provided for the use of other solvents, such as DMSO, water, or the culture medium. The detailed procedure for dissolving chemicals is shown in Appendix V. The X-VIVO™ 15 solutions described in paragraphs 18-23 are mixed 1:1 (v/v) with the cell suspensions prepared in a 96-well flat-bottom black plate (see paragraph 16).

18. The first test run is aimed to determine the cytotoxic concentration and to examine the skin sensitising potential of chemicals. Using X-VIVO™ 15, serial dilutions of the X-VIVO™ 15 stock solutions of the test chemicals are made at a dilution factor of two (see Appendix V) using a 96-well assay block (e.g. Costar Cat#EW-01729-03). Next, 50 µl/well of diluted solution is added to 50 µl of the cell suspension in a 96-well flat-bottom black plate. Thus for test chemicals that are soluble in X-VIVO™ 15, the final concentrations of the test chemicals range from 0.002 to 2 mg/mL (Appendix V). For test chemicals that are not soluble in X-VIVO™ 15 at 20 mg/mL, only dilution factors that range from 2 to 2¹⁰, are determined, although the actual final concentrations of the test chemicals remain uncertain and are dependent on the saturated concentration of the test chemicals in the X-VIVO™ 15 stock solution.

19. In subsequent test runs (i.e. the second, third, and fourth replicates), the X-VIVO™ 15 stock solution is made at the concentration 4 times higher than the concentration of cell viability 05 (CV05; the lowest concentration at which the Inh-GAPLA becomes <0.05) in the first experiment. If Inh-GAPLA does not decrease below 0.05 at the highest concentration in the first run, the X-VIVO™ 15 stock solution is made at the first run highest concentration. The concentration of CV05 is calculated by dividing the concentration of the stock solution in the first run by dilution factor for CV05 (X) (dilution factor CV05 (X)); the dilution factor required to dilute stock solution to CV05 (see Appendix V). For test substances not soluble in X-VIVO at 20 mg/ml, CV05 is determined by the concentration of the stock solution x 1/X. For run 2 to 4, a second stock solution is prepared as 4 x CV50 (Appendix V).

20. Serial dilutions of the X-VIVO™ 15 second stock solutions are made at a dilution factor of 1.5 using a 96-well assay block. Next, 50 µl/well of diluted solution is added to 50 µl of the cell suspension in the wells of a 96-well flat-bottom black plate. Each concentration of each test chemical should be tested in 4 wells. The samples are then mixed on a plate shaker and incubated for 16 hours at 37°C and 5% CO₂, after which the luciferase activity is measured as described below.

21. The solvent control is the mixture of 50 µL/well of X-VIVO™ 15 and 50 µL/well of cell suspension in RPMI-1640 containing 10% FBS.

22. The recommended positive control is 4-NBB. 20 mg of 4-NBB is prepared in a 1.5-mL microfuge tube, to which X-VIVO™ 15 is added up to 1 mL. The tube is vortexed vigorously and shaken on a rotor at a maximum speed of 8 rpm for at least 30 min. After centrifugation at 20,000g for 5 min, the supernatant is diluted by a factor of 4 with X-VIVO™ 15, and 500 µl of the diluted supernatant is transferred to a well in a 96-well assay block. The diluted supernatant is further diluted with X-VIVO™ 15 at factors of 2 and 4, and 50 µl of the solution is added to 50 µl of THP-G8 cell suspension in the wells of a 96-well flat-bottom black plate (Appendix VI). Each concentration of the positive control should be tested in 4 wells. The plate is agitated on a plateshaker, and incubated in a CO₂ incubator for 16 hours (37°C, 5% CO₂), after which the luciferase activity is measured as described in paragraph 29.

23. The recommended negative control is LA. 20 mg of LA prepared in a 1.5-mL microfuge tube, to which X-VIVO™ 15 is added up to 1 mL (20 mg/ mL). Twenty mg/mL of LA solution is diluted by a factor of 5 with X-VIVO™ 15 (4 mg/mL); 500 µl of this 4 mg/mL LA solution is transferred to a well of a 96-well assay block. This solution is diluted by a factor of 2 with X-VIVO™ 15 and then diluted again by a factor of 2 to produce 2 mg/mL and 1 mg/mL solutions. 50 µl of these 3 solutions and vehicle control (X-VIVO™ 15) are added to 50 µl of THP-G8 cell suspension in the wells of a 96-well flat-bottom black plate. Each concentration of the negative control is tested in 4 wells. The plate is agitated on a plateshaker and incubated in a CO₂ incubator for 16 hours (37°C, 5% CO₂), after which the luciferase activity is measured as described in paragraph 29.

24. Other suitable positive or negative controls may be used if historical data are available to derive comparable run acceptance criteria.
25. Care should be taken to avoid evaporation of volatile test chemicals and cross-contamination between wells by test chemicals, e.g. by sealing the plate prior to the incubation with the test chemicals.
26. The test chemicals and solvent control require 2 to 4 runs to derive a positive or negative prediction (see Table 2). Each run is performed on a different day with fresh X-VIVO™ 15 stock solution of test chemicals and independently harvested cells. Cells may come from the same passage.

Luciferase activity measurements

27. Luminescence is measured using a 96-well microplate luminometer equipped with optical filters, e.g. Phelios (ATTO, Tokyo, Japan), Tristan 941 (Berthold, Bad Wildbad, Germany) and the ARVO series (PerkinElmer, Waltham, MA, USA). The luminometer must be calibrated for each test to ensure reproducibility (19). Recombinant orange and red emitting luciferases are available for this calibration.
28. 100µL of pre-warmed Tripluc® Luciferase assay reagent (Tripluc) is transferred to each well of the plate containing the cell suspension treated with or without chemical. The plate is shaken for 10 min at an ambient temperature of about 20°C. The plate is placed in the luminometer to measure the luciferase activity. Bioluminescence is measured for 3 sec each in the absence (F0) and presence (F1) of the optical filter. Justification should be provided for the use of alternative settings, e.g. depending on the model of luminometer used.
29. Parameters for each concentration are calculated from the measured values, e.g. IL8LA, GAPLA, nIL8LA, Ind-IL8LA, Inh-GAPLA, the mean ±SD of IL8LA, the mean ±SD of GAPLA, the mean ±SD of nIL8LA, the mean ±SD of Ind-IL8LA, the mean ±SD of Inh-GAPLA, and the 95% confidence interval of Ind-IL8LA. Definitions of the parameters used in this paragraph are provided in Appendices I and IV, respectively.
30. Prior to measurement, colour discrimination in multi-colour reporter assays is generally achieved using detectors (luminometer and plate reader) equipped with optical filters, such as sharp-cut (long-pass or short-pass) filters or band-pass filters. The transmission coefficients of the filters for each bioluminescence signal colour should be calibrated prior to testing, per Appendix II.

DATA AND REPORTING

Data evaluation

31. Criteria for a positive/negative decision require that in each run:

- an IL-8 Luc assay prediction is judged positive if a test chemical has a Ind-IL8LA ≥ 1.4 and the lower limit of the 95% confidence interval of Ind-IL8LA ≥ 1.0
- an IL-8 Luc assay prediction is judged negative if a test chemical has a Ind-IL8LA < 1.4 and/or the lower limit of the 95% confidence interval of Ind-IL8LA < 1.0


Prediction model

32. Test chemicals that provide two positive results from among the 1st, 2nd, 3rd or 4th runs are identified as positives whereas those that give three negative results from among the 1st, 2nd, 3rd or 4th runs are identified as supposed negative (Table 2). Among supposed negative chemicals, chemicals that are dissolved at 20 mg/ml of X-VOVO™ 15 are judged as negative, while chemicals that are not dissolved at 20 mg/ml of X-VOVO™ 15 should not be considered (Figure 1).

Table 2. Criteria for identifying positive and supposed negative

1st run	2nd run	3rd run	4th run	Final prediction	
Positive	Positive	-	-	Positive	
	Negative	Positive	-	Positive	
		Negative	Positive	Negative	Supposed negative
	Negative	Positive	Positive	-	Positive
Negative			Positive	Positive	
Negative		Positive	Positive	Negative	Supposed negative
			Negative	Positive	Positive
		Negative	Positive	Negative	Supposed negative
			Negative	-	Supposed negative

Figure 1. Prediction model for final judgment

Acceptance criteria

33. The following acceptance criteria should be met when using the IL-8 Luc assay:

- Ind-IL8LA should be more than 5.0 at least in one concentration of the positive control, 4-NBB, in each run.
- Ind-IL8LA should be less than 1.4 at any concentration of the negative control, lactic acid, in each run.
- Data from plates for which the GAPLA of control wells with cells and Tripluc but without chemicals is less than 5 times of that of well containing test medium only (50 μ L/well of RPMI-1640 containing 10% FBS and 50 μ L/well of X-VIVO™ 15) should be rejected.
- Data from plates for which the Inh-GAPLA of all concentrations of the test or control chemicals is less than 0.05 should be rejected. In this case, the first test should be repeated so the highest final concentration of the repeated test is the lowest final concentration of the previous test.

TEST REPORT

34. The test report should include the following information:

Test chemicals

- Mono-constituent substance:
 - Chemical identification, such as IUPAC or CAS name(s), CAS number(s), SMILES or InChI code, structural formula, and/or other identifiers;
 - Physical appearance, water solubility, molecular weight, and additional relevant physicochemical properties, to the extent available;
 - Purity, chemical identity of impurities as appropriate and practically feasible, etc;

- Treatment prior to testing, if applicable (e.g. warming, grinding);
 - Solubility in X-VIVO™ 15. For chemicals that are insoluble in X-VIVO™ 15, whether precipitation or flotation are observed after centrifugation;
 - Concentration(s) tested;
 - Storage conditions and stability to the extent available;
 - Justification for choice of solvent/vehicle for each test chemical if X-VIVO™ 15 has not been used.
- Multi-constituent substance, UVCB and mixture:
- Characterisation as far as possible by e.g. chemical identity (see above), purity, quantitative occurrence and relevant physicochemical properties (see above) of the constituents, to the extent available;
 - Physical appearance, water solubility, and additional relevant physicochemical properties, to the extent available;
 - Molecular weight or apparent molecular weight in case of mixtures/polymers of known compositions or other information relevant for the conduct of the study;
 - Treatment prior to testing, if applicable (e.g. warming, grinding);
 - Solubility in X-VIVO™ 15. For chemicals that are insoluble in X-VIVO™ 15, whether precipitation or flotation are observed after centrifugation;
 - Concentration(s) tested;
 - Storage conditions and stability to the extent available.
 - Justification for choice of solvent/vehicle for each test chemical, if X-VIVO™ 15 has not been used.

Controls

- Positive control:
- Chemical identification, such as IUPAC or CAS name(s), CAS number(s), SMILES or InChI code, structural formula, and/or other identifiers;
 - Physical appearance, water solubility, molecular weight, and additional relevant physicochemical properties, to the extent available and where applicable;
 - Purity, chemical identity of impurities as appropriate and practically feasible, etc;
 - Treatment prior to testing, if applicable (e.g. warming, grinding);
 - Concentration(s) tested;
 - Storage conditions and stability to the extent available;
 - Reference to historical positive control results demonstrating suitable acceptance criteria, if applicable.

- Negative control:
 - Chemical identification, such as IUPAC or CAS name(s), CAS number(s), and/or other identifiers;
 - Purity, chemical identity of impurities as appropriate and practically feasible, etc;
 - Physical appearance, molecular weight, and additional relevant physicochemical properties in the case other negative controls than those mentioned in the Test Guideline are used and to the extent available;
 - Storage conditions and stability to the extent available;
 - Justification for choice of solvent for each test chemical.

Test method conditions

- Name and address of the sponsor, test facility and study director;
- Description of test method used;
- Cell line used, its storage conditions, and source (e.g. the facility from which it was obtained);
- Lot number and origin of FBC, supplier name, lot number of 96-well flat-bottom black plate, and lot number of Tripluc reagent;
- Passage number and cell density used for testing;
- Cell counting method used for seeding prior to testing and measures taken to ensure homogeneous cell number distribution;
- Luminometer used (e.g. model), including instrument settings, luciferase substrate used, and demonstration of appropriate luminescence measurements based on the control test described in Appendix II;
- The procedure used to demonstrate proficiency of the laboratory in performing the test method (e.g. by testing of proficiency substances) or to demonstrate reproducible performance of the test method over time.

Test procedure

- Number of replicates and runs performed;
- Test chemical concentrations, application procedure and exposure time (if different from those recommended);
- Description of evaluation and decision criteria used;
- Description of study acceptance criteria used;
- Description of any modifications of the test procedure.

Results

- Measurements of IL8LA and GAPLA;
- Calculations for nIL8LA, Ind-IL8LA, and Inh-GAPLA;
- The 95% confidence interval of Ind-IL8LA;
- A graph depicting dose-response curves for induction of luciferase activity and viability;
- Description of any other relevant observations, if applicable.

Discussion of the results

- Discussion of the results obtained with the IL-8 Luc assay;
- Consideration of the assay results in the context of an IATA, if other relevant information is available.

Conclusion

LITERATURE

1. Takahashi T, Kimura Y, Saito R, Nakajima Y, Ohmiya Y, Yamasaki K, and Aiba S. 2011. An in vitro test to screen skin sensitizers using a stable THP-1-derived IL-8 reporter cell line, THP-G8. *Toxicol Sci* 124:359-69.
2. OECD (2017), Validation report for the international validation study on the IL-8 Luc assay as a test evaluating the skin sensitizing potential of chemicals conducted by the IL-8 Luc Assay. Series on Testing and Assessment No. 267, [ENV/JM/MONO\(2017\)19](#). Organisation for Economic Cooperation and Development, Paris. Available at: [<http://www.oecd.org/env/ehs/testing/series-testing-assessment-publications-number.htm>].
3. OECD (2017), Report of the Peer Review Panel for the IL-8 Luciferase (IL-8 Luc) Assay for in vitro skin sensitisation. Series on Testing and Assessment No. 258, [ENV/JM/MONO\(2017\)20](#). Organisation for Economic Cooperation and Development, Paris. Available at: [<http://www.oecd.org/env/ehs/testing/series-testing-assessment-publications-number.htm>].
2. OECD (2016), Guidance Document On The Reporting Of Defined Approaches And Individual Information Sources To Be Used Within Integrated Approaches To Testing And Assessment (IATA) For Skin Sensitisation, Series on Testing & Assessment No. 256, [ENV/JM/MONO\(2016\)29](#). Organisation for Economic Cooperation and Development, Paris. Available at: [<http://www.oecd.org/env/ehs/testing/series-testing-assessment-publications-number.htm>].
3. van der Veen JW, Rorije E, Emter R, Natsch A, van Loveren H, and Ezendam J. 2014. Evaluating the performance of integrated approaches for hazard identification of skin sensitizing chemicals. *Regul Toxicol Pharmacol* 69:371-9.
4. Kimura Y, Fujimura C, Ito Y, Takahashi T, Nakajima Y, Ohmiya Y, and Aiba S. 2015. Optimization of the IL-8 Luc assay as an in vitro test for skin sensitization. *Toxicol In Vitro* 29:1816-30.
5. Urbisch D, Mehling A, Guth K, Ramirez T, Honarvar N, Kolle S, Landsiedel R, Jaworska J, Kern PS, Gerberick F, et al. 2015. Assessing skin sensitization hazard in mice and men using non-animal test methods. *Regul Toxicol Pharmacol* 71:337-51.
6. Ashikaga T, Sakaguchi H, Sono S, Kosaka N, Ishikawa M, Nukada Y, Miyazawa M, Ito Y, Nishiyama N, and Itagaki H. 2010. A comparative evaluation of in vitro skin sensitisation tests: the human cell-line activation test (h-CLAT) versus the local lymph node assay (LLNA). *Alternatives to laboratory animals* : *ATLA* 38:275-84.
7. Patlewicz G, Casati S, Basketter DA, Asturiol D, Roberts DW, Lepoittevin J-P, Worth A and Aschberger K (2016) Can currently available non-animal methods detect pre and pro haptens relevant for skin sensitisation? *Regul Toxicol Pharmacol*, online.
8. Thorne N, Inglese J, and Auld DS. 2010. Illuminating insights into firefly luciferase and other bioluminescent reporters used in chemical biology. *Chem Biol* 17:646-57.
9. OECD (2016), *Test No. 455: Performance-Based Test Guideline for Stably Transfected Transactivation In Vitro Assays to Detect Estrogen Receptor Agonists and Antagonists*, OECD Publishing, Paris. <http://dx.doi.org/10.1787/9789264265295-en>.

10. Viviani V, Uchida A, Suenaga N, Ryufuku M, and Ohmiya Y. 2001. Thr226 is a key residue for bioluminescence spectra determination in beetle luciferases. *Biochem Biophys Res Commun* 280:1286-91.
11. Viviani VR, Bechara EJ, and Ohmiya Y. 1999. Cloning, sequence analysis, and expression of active *Phrixothrix* railroad-worms luciferases: relationship between bioluminescence spectra and primary structures. *Biochemistry* 38:8271-9.
12. Nakajima Y, Kimura T, Sugata K, Enomoto T, Asakawa A, Kubota H, Ikeda M, and Ohmiya Y. 2005. Multicolor luciferase assay system: one-step monitoring of multiple gene expressions with a single substrate. *Biotechniques* 38:891-4.
13. OECD (2017), To be published - Performance Standards for the assessment of proposed similar or modified *in vitro* skin sensitisation IL-8 luc test methods. OECD Environment, Health and Safety Publications, Series on Testing and Assessment. OECD, Paris, France.
14. OECD (2005), Guidance Document the Validation and International Acceptance of New or Updated Test Methods for Hazard Assessment. OECD Environment, Health and Safety publications, OECD Series on Testing and Assessment No.34. OECD, Paris, France.
15. OECD (2017), Draft Guidance document: Good *In Vitro* Method Practices (GIVIMP) for the Development and Implementation of *In Vitro* Methods for Regulatory Use in Human Safety Assessment. Organisation for Economic Cooperation and Development, Paris. Available at: [http://www.oecd.org/env/ehs/testing/OECD%20Draft%20GIVIMP_v05%20-%20clean.pdf].
16. JaCVAM (2016), IL-8 Luc assay protocol, Available at. http://www.jacvam.jp/en_effort/effort02.html.
17. Niwa K, Ichino Y, Kumata S, Nakajima Y, Hiraishi Y, Kato D, Viviani VR, and Ohmiya Y. 2010. Quantum yields and kinetics of the firefly bioluminescence reaction of beetle luciferases. *Photochem Photobiol* 86:1046-9.
18. OECD (2012), The Adverse Outcome Pathway for Skin Sensitisation Initiated by Covalent Binding to Proteins, Part 1: Scientific Evidence. OECD Environment, Health and Safety Publications, Series on Testing and Assessment No. 168. OECD, Paris, France.
19. United Nations (2015), Globally Harmonized System of Classification and Labelling of Chemicals (GHS). Sixth revised edition. New York & Geneva: United Nations Publications. ISBN: 978-92-1-117006-1. Available at: http://www.unece.org/trans/danger/publi/ghs/ghs_rev05/05files_e.html.

APPENDIX I

DEFINITIONS

Accuracy: The closeness of agreement between test method results and accepted reference values. It is a measure of test method performance and one aspect of relevance. The term is often used interchangeably with concordance to mean the proportion of correct outcomes of a test method (16).

AOP (Adverse Outcome Pathway): Sequence of events from the chemical structure of a target chemical or group of similar chemicals through the molecular initiating event to an *in vivo* outcome of interest (20).

CV05: Cell viability 05. Minimum concentration at which chemicals show less than 0.05 of Inh-GAPLA.

FInSLO-LA: Abbreviation used in the validation report and in previous publications regarding the IL-8 Luc assay to refer to Ind-IL8LA. See Ind-IL8LA for definition.

GAPLA: Luciferase Activity of Stable Luciferase Red (SLR) ($\lambda_{\max} = 630$ nm), regulated by GAPDH promoter and demonstrates cell viability and viable cell number.

Hazard: Inherent property of an agent or situation having the potential to cause adverse effects when an organism, system or (sub) population is exposed to that agent.

IATA (Integrated Approach to Testing and Assessment): A structured approach used for hazard identification (potential), hazard characterisation (potency) and/or safety assessment (potential/potency and exposure) of a chemical or group of chemicals, which strategically integrates and weights all relevant data to inform regulatory decision regarding potential hazard and/or risk and/or the need for further targeted and therefore minimal testing.

II-SLR-LA: Abbreviation used in the validation report and in previous publications regarding the IL-8 Luc assay to refer to Inh-GAPLA. See Inh-GAPLA for definition

IL-8 (Interleukin-8): A cytokine derived from endothelial cells, fibroblasts, keratinocytes, macrophages, and monocytes that causes chemotaxis of neutrophils and T-cell lymphocytes.

IL8LA: Luciferase Activity of Stable Luciferase Orange (SLO) ($\lambda_{\max} = 580$ nm), regulated by IL-8 promoter.

Ind-IL8LA: Fold induction of IL8LA. It is obtained by dividing the nIL8LA of THP-G8 cells treated with chemicals by that of non-stimulated THP-G8 cells and represents the induction of IL-8 promoter activity by chemicals.

Inh-GAPLA: Inhibition of GAPLA. It is obtained by dividing GAPLA of THP-G8 treated with chemicals with GAPLA of non-treated THP-G8 and represents cytotoxicity of chemicals.

Minimum induction threshold (MIT): the lowest concentration at which a chemical satisfies the positive criteria

Mixture: A mixture or a solution composed of two or more substances in which they do not react.

Mono-constituent substance: A substance, defined by its quantitative composition, in which one main constituent is present to at least 80% (w/w).

Multi-constituent substance: A substance, defined by its quantitative composition, in which more than one of the main constituents is present in a concentration $\geq 10\%$ (w/w) and $< 80\%$ (w/w). A multi-constituent substance is the result of a manufacturing process. The difference between mixture and multi-constituent substance is that a mixture is obtained by blending of two or more substances without chemical reaction. A multi-constituent substance is the result of a chemical reaction.

nIL8LA: The SLO luciferase activity reflecting IL-8 promoter activity (IL8LA) normalised by the SLR luciferase activity reflecting GAPDH promoter activity (GALPA). It represents IL-8 promoter activity after considering cell viability or cell number.

nSLO-LA: Abbreviation used in the validation report and in previous publications regarding the IL-8 Luc assay to refer to nIL8LA. See nIL8LA for definition

Positive control: A replicate containing all components of a test system and treated with a substance known to induce a positive response. To ensure that variability in the positive control response across time can be assessed, the magnitude of the positive response should not be excessive.

Pre-haptens: Chemicals which become sensitisers through abiotic transformation.

Pro-haptens: Chemicals requiring enzymatic activation to exert skin sensitisation potential.

Relevance: Description of relationship of the test to the effect of interest and whether it is meaningful and useful for a particular purpose. It is the extent to which the test correctly measures or predicts the biological effect of interest. Relevance incorporates consideration of the accuracy (concordance) of a test method (16).

Reliability: Measures of the extent that a test method can be performed reproducibly within and between laboratories over time, when performed using the same protocol. It is assessed by calculating intra- and inter-laboratory reproducibility and intra-laboratory repeatability (16).

Run: A run consists of one or more test chemicals tested concurrently with a solvent/vehicle control and with a positive control.

Sensitivity: The proportion of all positive/active chemicals that are correctly classified by the test. It is a measure of accuracy for a test method that produces categorical results, and is an important consideration in assessing the relevance of a test method (16).

SLO-LA: Abbreviation used in the validation report and in previous publications regarding the IL-8 Luc assay to refer to IL8LA. See IL8LA for definition.

SLR-LA: Abbreviation used in the validation report and in previous publications regarding the IL-8 Luc assay to refer to GAPLA. See GAPLA for definition.

Solvent/vehicle control: An untreated sample containing all components of a test system except of the test chemical, but including the solvent/vehicle that is used. It is used to establish the baseline response

for the samples treated with the test chemical dissolved or stably dispersed in the same solvent/vehicle. When tested with a concurrent medium control, this sample also demonstrates whether the solvent/vehicle interacts with the test system.

Specificity: The proportion of all negative/inactive chemicals that are correctly classified by the test. It is a measure of accuracy for a test method that produces categorical results and is an important consideration in assessing the relevance of a test method (16).

Substance: Chemical elements and their compounds in the natural state or obtained by any production process, including any additive necessary to preserve the stability of the product and any impurities deriving from the process used, but excluding any solvent which may be separated without affecting the stability of the substance or changing its composition.

Surfactant: Also called surface-active agent, this is a substance, such as a detergent, that can reduce the surface tension of a liquid and thus allow it to foam or penetrate solids; it is also known as a wetting agent. (TG437)

Test chemical: The term "test chemical" is used to refer to what is being tested.

THP-G8: An IL-8 reporter cell line used in IL-8 Luc assay. The human macrophage-like cell line THP-1 was transfected the SLO and SLR luciferase genes under the control of the IL-8 and GAPDH promoters, respectively.

United Nations Globally Harmonized System of Classification and Labeling of Chemicals (UN GHS): A system proposing the classification of chemicals (substances and mixtures) according to standardised types and levels of physical, health and environmental hazards, and addressing corresponding communication elements, such as pictograms, signal words, hazard statements, precautionary statements and safety data sheets, so that to convey information on their adverse effects with a view to protect people (including employers, workers, transporters, consumers and emergency responders) and the environment (21).

UVCB: substances of unknown or variable composition, complex reaction products or biological materials.

Valid test method: A test method considered to have sufficient relevance and reliability for a specific purpose and which is based on scientifically sound principles. A test method is never valid in an absolute sense, but only in relation to a defined purpose.

APPENDIX II

**PRINCIPLE OF MEASUREMENT OF LUCIFERASE ACTIVITY AND DETERMINATION
OF THE TRANSMISSION COEFFICIENTS OF OPTICAL FILTER FOR SLO AND SLR**

MultiReporter Assay System -Tripluc- can be used with a microplate-type luminometer with a multi-colour detection system, which can equip an optical filter (e.g. Phelios AB-2350 (ATTO), ARVO (PerkinElmer), Tristar LB941 (Berthold)). The optical filter used in measurement is 600–620 nm long or short pass filter, or 600–700 nm band pass filter.

(1) Measurement of two-colour luciferases with an optical filter.

This is an example using Phelios AB-2350 (ATTO). This luminometer is equipped with a 600 nm long pass filter (R60 HOYA Co.), 600 nm LP, Filter 1) for splitting SLO ($\lambda_{\max} = 580$ nm) and SLR ($\lambda_{\max} = 630$ nm) luminescence.

To determine transmission coefficients of the 600 nm LP, first, using purified SLO and SLR luciferase enzymes, measure i) the intensity of SLO and SLR bioluminescence intensity without filter (F0), ii) the SLO and SLR bioluminescence intensity that passed through 600 nm LP (Filter 1), and iii) calculate the transmission coefficients of 600 nm LP for SLO and SLR listed below.

Transmission coefficients	Abbreviation	Definition
SLO Filter 1 Transmission coefficients	$\kappa_{O_{R60}}$	The filter's transmission coefficient for the SLO
SLR Filter 1 Transmission coefficients	$\kappa_{R_{R60}}$	The filter's transmission coefficient for the SLR

When the intensity of SLO and SLR in test sample are defined as O and R, respectively, i) the intensity of light without filter (all optical) F0 and ii) the intensity of light that transmits through 600 nm LP (Filter 1) F1 are described as below.

$$F0 = O + R$$

$$F1 = \kappa_{O_{R60}} \times O + \kappa_{R_{R60}} \times R$$

These formulas can be rephrased as follows:

$$\begin{pmatrix} F0 \\ F1 \end{pmatrix} = \begin{pmatrix} 1 & 1 \\ \kappa_{O_{R60}} & \kappa_{R_{R60}} \end{pmatrix} \begin{pmatrix} O \\ R \end{pmatrix}$$

Then using calculated transmittance factors ($\kappa_{O_{R60}}$ and $\kappa_{R_{R60}}$) and measured F0 and F1, you can calculate O and R-value as follows:

$$\begin{pmatrix} O \\ R \end{pmatrix} = \begin{pmatrix} 1 & 1 \\ \kappa_{O_{R60}} & \kappa_{R_{R60}} \end{pmatrix}^{-1} \begin{pmatrix} F0 \\ F1 \end{pmatrix}$$

Materials and methods for determining transmittance factor

(1) Reagents

· Single purified luciferase enzymes:

Lyophilised purified SLO enzyme

Lyophilised purified SLR enzyme

(which for the validation work were obtained from GPC Lab. Co. Ltd., Tottori, Japan with THP-G8 cell line)

· Assay reagent:

Tripluc[®] Luciferase assay reagent (for example from TOYOBO Cat#MRA-301)

· Medium: for luciferase assay (30 ml, stored at 2 – 8°C)

Reagent	Conc.	Final conc. in medium	Required amount
RPMI-1640	-	-	27 ml
FBS	-	10 %	3 ml

(2) Preparation of enzyme solution

Dissolve lyophilised purified luciferase enzyme in tube by adding 200 µl of 10 ~ 100 mM Tris/HCl or Hepes/HCl (pH 7.5 ~ 8.0) supplemented with 10% (w/v) glycerol, divide the enzyme solution into 10 µl aliquots in 1.5 ml disposable tubes and store them in a freezer at -80°C. The frozen enzyme solution can be used for up to 6 months. When used, add 1 ml of medium for luciferase assay (RPMI-1640 with 10% FBS) to each tube containing the enzyme solutions (diluted enzyme solution) and keep them on ice to prevent deactivation.

(3) Bioluminescence measurement

Thaw Tripluc[®] Luciferase assay reagent (Tripluc) and keep it at room temperature either in a water bath or at ambient air temperature. Power on the luminometer 30 min before starting the measurement to allow the photomultiplier to stabilise. Transfer 100 µl of the diluted enzyme solution to a black 96 well plate (flat bottom) (the SLO reference sample to #B1, #B2, #B3, the SLR reference sample to #D1, #D2, #D3). Then, transfer 100 µl of pre-warmed Tripluc to each well of the plate containing the diluted enzyme solution using a pipetman. Shake the plate for 10 min at room temperature (about 25°C) using a plate shaker. Remove bubbles from the solutions in wells if they appear. Place the plate in the luminometer to measure the luciferase activity. Bioluminescence is measured for 3 sec each in the absence (F0) and presence (F1) of the optical filter.

Transmission coefficient of the optical filter was calculated as follows:

Transmission coefficient (SLO ($\kappa_{O_{R60}}$))= (#B1 of F1+ #B2 of F1+ #B3 of F1) / (#B1 of F0+ #B2 of F0+ #B3 of F0)

Transmission coefficient (SLR ($\kappa_{R_{60}}$))= (#D1 of F1+ #D2 of F1+ #D3 of F1) / (#D1 of F0+ #D2 of F0+ #D3 of F0)

Calculated transmittance factors are used for all the measurements executed using the same luminometer.

Quality control of equipment

The procedures described in the IL-8 Luc protocol should be used (18).

APPENDIX III

PROFICIENCY SUBSTANCES

Prior to routine use of the test method described in this Annex to Test Guideline 442E, laboratories should demonstrate technical proficiency by obtaining the expected IL-8 Luc assay prediction for the 9 substances recommended in Table 1 and by obtaining values that fall within the respective reference range for at least 8 out of the 9 proficiency substances (selected to represent the range of responses for skin sensitisation hazards). Other selection criteria were that the substances are commercially available, and that high-quality *in vivo* reference data as well as high quality *in vitro* data generated with the IL-8 Luc assay are available. Also, published reference data are available for the IL-8 Luc assay (6) (1).

Table 1: Recommended substances for demonstrating technical proficiency with the IL-8 Luc assay

Proficiency substances	CAS no.	State	Solubility in X- VIVO15 at 20 mg/mL	<i>In vivo</i> prediction ¹	IL-8 Luc prediction ²	Reference range (µg/mL) ³	
						CV05 ⁴	IL-8 Luc MIT ⁵
2,4-Dinitrochlorobenzene	97-00-7	Solid	Insoluble	Sensitiser (Extreme)	Positive	2.3-3.9	0.5-2.3
Formaldehyde	50-00-0	Liquid	Soluble	Sensitiser (Strong)	Positive	9-30	4-9
2-Mercaptobenzothiazole	149-30-4	Solid	Insoluble	Sensitiser (Moderate)	Positive	250-290	60-250
Ethylenediamine	107-15-3	Liquid	Soluble	Sensitiser (Moderate)	Positive	500-700	0.1-0.4
Ethyleneglycol dimethacrylate	97-90-5	Liquid	Insoluble	Sensitiser (Weak)	Positive	>2000	0.04-0.1
4-Allylanisole (Estragol)	140-67-0	Liquid	Insoluble	Sensitiser (Weak)	Positive	>2000	0.01-0.07
Streptomycin sulphate	3810-74-0	Solid	Soluble	Non- sensitiser	Negative	>2000	>2000
Glycerol	56-81-5	Liquid	Soluble	Non- sensitiser	Negative	>2000	>2000
Isopropanol	67-63-0	Liquid	Soluble	Non- sensitiser	Negative	>2000	>2000

Abbreviations: CAS no. = Chemical Abstracts Service Registry Number

¹ The *in vivo* potency is derived using the criteria proposed by ECETOC (19).

² Based on historical observed values (1) (6).

³ CV05 and IL-8 Luc MIT were calculated using water solubility given by EPI Suite™.

⁴ CV05: the minimum concentration at which chemicals show less than 0.05 of Inh-GAPLA.

⁵ MIT: the lowest concentrations at which a chemical satisfies the positive criteria.

APPENDIX IV

INDEXES AND JUDGMENT CRITERIA

nIL8LA (nSLO-LA)

The j-th repetition (j = 1-4) of the i-th concentration (i = 0-11) is measured for IL8LA (SLO-LA) and GAPLA (SLR-LA) respectively. The normalised IL8LA, referred to as nIL8LA (nSLO-LA), and is defined as:

$$nIL8LA_{ij} = IL8LA_{ij} / GAPLA_{ij}$$

This is the basic unit of measurement in this assay.

Ind-IL8LA (FInSLO-LA)

The fold increase of the averaged nIL8LA (nSLO-LA) for the repetition on the i-th concentration compared with it at the 0 concentration, Ind-IL8LA, is the primary measure of this assay. This ratio is written by the following formula:

$$Ind-IL8LA_i = \left\{ (1/4) \times \sum_j nIL8LA_{ij} \right\} / \left\{ (1/4) \times \sum_j nIL8LA_{0j} \right\}$$

The lead laboratory has proposed that a value of 1.4 corresponds to a positive result for the tested chemical. This value is based on the investigation of the historical data of the lead laboratory. Data management team then used this value through all the phases of validation study. The primary outcome, Ind-IL8LA, is the ratio of 2 arithmetic means as shown in equation.

95% confidence interval (95% CI)

The 95% confidence interval (95% CI) based on the ratio can be estimated to show the precision of this primary outcome measure. The lower limit of the 95% CI ≥ 1 indicates that the nIL8LA with the i-th concentration is significantly greater than that with solvent control. There are several ways to construct the 95% CI. We used the method known as Fieller's theorem in this study. This 95% confidence interval theorem is obtained from the following formula:

$$\left[\frac{-B - \sqrt{B^2 - 4AC}}{2A}, \frac{-B + \sqrt{B^2 - 4AC}}{2A} \right],$$

where $A = \bar{x}_0^2 - t_{0.975(\nu)}^2 \times \frac{sd_0^2}{n_0}$, $B = -2 \times \bar{x} \times \bar{y}$, $C = \bar{y}_i^2 - t_{0.975(\nu)}^2 \times \frac{sd_{y_i}^2}{n_{y_i}}$, and $n_0 = 4$,

$$\bar{x}_0 = (1/n_0) \times \sum_j nIL8LA_{0j}, \quad sd_0^2 = \{1/(n_0 - 1)\} \times \sum_j (nIL8LA_{0j} - \bar{x}_0)^2,$$

$$n_{y_i} = 4, \quad \bar{y}_i = (1/n_{y_i}) \times \sum_j (nIL8LA_{ij}), \quad sd_{y_i}^2 = \{1/(n_{y_i} - 1)\} \times \sum_j (nIL8LA_{ij} - \bar{y}_i)^2.$$

$t_{0.975(\nu)}$ is 97.5 percentile of the central t distribution with the ν of the degree of freedom, where

$$v = \left(\frac{sd_0^2}{n_0} + \frac{sd_{y_i}^2}{n_{y_i}} \right)^2 / \left\{ \left(\frac{sd_0^2}{n_0} \right)^2 / (n_0 - 1) + \left(\frac{sd_{y_i}^2}{n_{y_i}} \right)^2 / (n_{y_i} - 1) \right\}.$$

Inh-GAPLA (II-SLR-LA)

The Inh-GAPLA is a ratio of the averaged GAPLA (SLR-LA) for the repetition of the i-th concentration compared with that with solvent control, and this is written by


$$\text{Inh-GAPLA}_i = \left\{ (1/4) \times \sum_j \text{GAPLA}_{ij} \right\} / \left\{ (1/4) \times \sum_j \text{GAPLA}_{0j} \right\}.$$

Since the GAPLA is the denominator of the nIL8LA, an extremely small value causes large variation in the nIL8LA. Therefore, Ind-IL8LA values with an extremely small value of Inh-GAPLA (less than 0.05) might be considered poor precision.



APPENDIX V

THE SCHEME OF THE METHODS TO DISSOLVE CHEMICALS FOR THE IL-8 LUC ASSAY.




(a) For chemicals dissolved in X-VIVO™ 15 at 20 mg/mL

Determine the highest concentration of the following experiments □

(b) For chemicals insoluble in X-VIVO™ 15 at 20 mg/mL

APPENDIX VI

THE SCHEME OF THE METHOD TO DISSOLVE 4-NBB FOR THE POSITIVE CONTROL OF THE IL-8 LUC ASSAY.

The positive control : 4-NBB (insoluble in X-VIVO™ 15) □