

1
2
3
4
5
6
7
8
9
10
11
12
13
14
15
16
17
18
19
20
21
22
23
24
25
26
27
28
29
30
31
32
33
34
35
36
37
38
39
40
41
42

Hand1-Luc Embryonic Stem Cell Test (Hand1-Luc EST)

Draft protocol ver. 09 E

March, 9th, 2016

Environmental Health Science Laboratory,
Sumitomo Chemical Co., Ltd.

Noriyuki Suzuki, Ph.D.
Florian Le Coz, M.Sc.
Hirohisa Nagahori, Ph.D.
Koichi Saito, Ph.D.

TABLE OF CONTENTS

1		
2		
3	INTRODUCTION	3
4		
5	MATERIALS AND METHODS	3
6	EXPERIMENT 1	3
7	1-1 Measurement of actual value of pKa or <i>in silico</i> method	3
8		
9	EXPERIMENT 2	4
10	2-1 Measurement of ID₅₀ and IC₅₀ values in the Hand1-Luc EST	4
11	A. Materials	4
12	<i>Cells (Hand1-ES cells).....</i>	4
13	<i>Reagents and preparation methods.....</i>	4
14	<i>Equipments.....</i>	6
15	<i>Expendable supplies.....</i>	6
16	B. Experimental Methods	7
17	<i>Thawing of Hand1-ES cells.....</i>	7
18	<i>Maintenance of cells</i>	7
19	<i>Preparation of test chemicals</i>	8
20	<i>Serial dilution of test chemicals and arrangement of plate format.....</i>	9
21	<i>Dose finding study</i>	10
22	<i>Definitive study</i>	13
23	<i>Data analysis</i>	13
24	<i>Acceptance criteria</i>	14
25		
26	EXPERIMENT 3 (optional prediction with PM2)	15
27	3-1 Measurement of metabolic stability of test chemicals	15
28		
29	JUDGMENT (Prediction of test chemicals)	16
30		
31	OPTION.....	17
32		
33	UPDATE RECORD	19
34		

INTRODUCTION

Hand1-Luc Embryonic Stem Cell Test (Hand1-Luc EST) is a novel short term test for predicting embryonic toxic chemicals using transgenic engineering mouse embryonic stem (ES) cells with consideration of pharmacokinetic property of test chemicals. This protocol describes how to prepare, determine the pharmacokinetic property and analyze cytotoxicity and differentiation toxicity of chemicals.

MATERIALS AND METHODS

~~EXPERIMENT 1~~

~~1-1 Measurement of actual value of pKa or in silico methods~~

~~Actual value can be measured by titration, spectrophotometric and conductometric methods. Detailed protocol is shown in OECD Test Guideline (112,~~

~~<http://www.oecd.org/env/ehs/testing/oecdguidelinesforthetestingofchemicals.htm>~~

~~The value can be calculated by commercial software or available tools on the internet such as ADMET Predictor, ACD/labs or Mervin (<http://www.chemaxon.com/marvin/sketch/index.jsp>).~~

EXPERIMENT 2

2-1 Measurement of ID₅₀ and IC₅₀ values in the Hand1-Luc EST

Schema of the Hand1-Luc EST procedure is described in Figure 1.




Figure.1 Outline of the Hand1-Luc EST

A. Materials

Cells (Hand1-ES)

Transgenic mouse embryonic stem (ES) cell, named as Hand1-ES, containing firefly luciferase genes under the control of Hand1 (heart and neural crest derivatives expressed transcript 1) promoters were established by Sumitomo Chemical Co., Ltd. as a previous method with minor modification (Suzuki et al., Toxicol. Sci., 124, 460–471, 2011). An embryonic stem cell line (named as KOB1-ES) from C57/B6 mouse established by Sumitomo Chemical Co. Ltd. was used for development of the Hand1-ES cells.

Reagents and preparation methods

For maintenance of Hand1-ES cells

- StemMedium (Serum Free Media for Mouse ES Cell) (DS Pharma Biomedical, Cat. DSRK100)

▲Note **Dissolved Dispensed** medium is stored at -20 °C (in an appropriate volume in tubes).

- Gelatin from porcine, type A (SIGMA, Cat. G2500)
- ESGRO mLIF (10⁷ units/mL) (Millipore, Cat. ESG1107)

- 1 – Distilled water (SIGMA, Cat. W3500-500mL)
 2 – 0.25 % Trypsin /1mM EDTA solution (Nacalai tesque, Cat.35554-64)
 3 – 50 mg/mL Geneticin Disulfate (G418) solution (Nacalai tesque, Cat.
 4 16513-84)
 5 – Trypan blue solution (SIGMA, Cat.T8154-20ML or equivalents, regardless
 6 of makers, are acceptable)

7
 8 For chemical exposure, positive control, solvents

- 9 – 5-Fluorouracil (5-FU) (SIGMA, Cat. F6627-1G, purity 99 %)
 10 – Dimethyl sulfoxide (DMSO) (MP Biomedicals, Inc. (WAKO Pure Chemicals),
 11 Cat. 594-03771)
 12 – PBS(-) (Invitrogen, Cat No. 10010-023 (500mL))

13
 14 For differentiation or cytotoxicity assay

- 15 – DMEM (Invitrogen, Cat. 11965-092)
 16 – FBS (Hyclone, Cat 533-90935, Lot. KTC30729)
 17 **▲Note** The quality of the serum is very important for mESC growth
 18 and differentiation. Several lots of serum have to be tested before
 19 performing the assay.
 20 – GlutaMaxI solution (100xconc.) (Invitrogen, Cat. 35050-061)
 21 – 100 mM Non-Essential Amino Acid solution (NEAA) (Invitrogen, Cat.
 22 11140-050)
 23 – Penicillin-Streptomycin, liquid (PS solution) (Nacalai tesque, Cat.
 24 26252-94)
 25 – 2-Mercaptoethanol (2-ME) (Nacalai tesque, Cat. 21438-82)
 26 – Steady-Glo® Luciferase Assay System (Promega, Cat No. E2510, E2520
 27 and E2550)
 28 **▲Note** Dissolved reagents stored at -80 °C (in an appropriate volume in
 29 tubes).
 30 – CellTiter-Fluor™ Cell Viability Assay (Promega, Cat No. G6080, G6081 and
 31 G6082)

32
 33 Maintenance medium (Culture medium for maintenance of Hand1-ES cells)

Content	Volume
StemMedium	10 mL
100mM 2-ME solution	10 μL
ESGRO (10 ⁷ unit/mL)	1 μL
50 mg/mL G418 solution	20 μL
Note	
- Supplemented medium is stored for no longer than 1 week at 4 °C.	
- Before use, supplemented medium should be prewarmed to 37 °C in a water bath or incubator.	

34
 35 Assay medium

Content	Volume
---------	--------

DMEM	82 mL
FBS (Heat inactivated)	15 mL
GlutaMaxI solution (100x conc.)	1 mL
100mM NEAA solution	1 mL
PS solution	1 mL
100mM 2-ME solution	100 μ L
Note	
- Supplemented medium is stored for no longer than 2 weeks at 4 °C.	
- Before use, supplemented medium should be prewarmed to 37 °C in a water bath or incubator.	

1

2

0.1% Gelatin solution

Content	Volume
Gelatin	0.5 g
Distilled water	500 mL
Note	
- Autoclaved solution is stored for no longer than 2 months at 4 °C.	
- 60mm dishes should be coated with gelatin solution (37°C, over 30min).	

3

4

2-Mercaptoethanol (2-ME) solution

Content	Volume
2-ME (14.3 M solution)	70 μ L
Distilled water	10 mL
Note	
Stock solutions are stored at -80 °C or prepared at time of use.	

5

6

Fetal Bovine Serum, Heat-Inactivated

Content
FBS
Note
Incubate at 56 °C for 30 min and store at -20 °C (appropriate volume in tubes).
Dissolved FBS is kept at 4 °C and use within 1 month.

7

8

9

Equipments

- 10 – Measuring device of micro-plate type luminometer
- 11 – Measuring device of micro-plate type fluorometer
- 12 – Plate shaker
- 13 – Water bath
- 14 – CO₂ incubator (5%, 37°C)
- 15 – Single and 8 channel micropipette
- 16 – Appropriate cell counter (hemocytometer)

17

18

Expendable supplies

19

20

- 60 mm Cell culture dish (BD Falcon, Cat.353004)
- Storage Plate 96 well Round Bottom (Corning, Cat.3359)

- 1 – 96 well Assay Blocks 2mL/well (Corning, Cat.3960)
- 2 – PrimeSurface® 96well white plate (Sumitomo Bake Lite, Co., Ltd.
3 MS-9096W)
- 4 – Plate seal (for PrimeSurface® 96well white plate) (Watson, Cat.
5 547-KTS-HC)
- 6 – Plate seal “TopSeal-A₊” (PerkinElmer No. 6050185)
- 7 – Cell strainer (40µm, BD Falcon, Cat.352340)
- 8 – Reservoir
- 9 – Pipette

10

11 **B. Experimental Methods**

12

13 *Thawing of Hand1-ES cells*

14

- 15 1. Coat the 60 mm dishes with the 0.1 % gelatin solution. Add 5 ml of the
16 gelatin solution in the dishes and incubate at 37 °C during 30 min.
- 17 2. Thaw frozen cells in 37 °C water bath, and add to the 9 mL of
18 maintenance medium.
- 19 3. Centrifuge the tube at 900 to 1,400 rpm for 5 min at room temperature,
20 discard supernatant and resuspend the cells in 5 mL of maintenance
21 medium.
- 22 4. Count cell number using aliquot of solution.
- 23 5. Remove the gelatin solution from the dishes (step 1) and seed the
24 cells at a concentration of approx. 0.5 to 1.0×10^6 cells/ gelatin-coated
25 60 mm dish in 5mL of maintenance medium.
- 26 6. Incubate the cells at 37°C in a humidified atmosphere with 5% CO₂.

27

28 *Maintenance of cells*

29

30 Grown cells should be passaged 2 or 3 days after thawing.

- 31 1. Prepare the 60 mm dishes coated with the 0.1 % gelatin solution
32 before maintenance of cells.
- 33 2. Remove the supernatant from culture dish, wash with 5 mL of PBS(-).
- 34 3. Add 2 mL of 0.25% trypsin/1 mM EDTA and remove immediately.
35 Then, cells are incubated at 37 °C for 1 to 2 min.
- 36 4. Add to 2 mL of maintenance medium to the dish, suspend the cells
37 using micropipette (1000 µL).
- 38 5. Count cell number
- 39 6. Remove the gelatin solution from the dishes (step 1)
- 40 7. For assay, seed the cells at a concentration of 2.0×10^6 cells/
41 gelatin-coated 60 mm dish in 5mL of maintenance medium. Incubate
42 the cells at 37°C with 5% CO₂ for 1 day.
- 43 8. For passage, seed the cells at a concentration of 0.2, 0.5 and
44 1.0×10^6 cells/ gelatin-coated 60 mm dish in 5mL of maintenance
45 medium. Incubate the cells at 37°C with 5% CO₂ for 2 or 3 days. Cells
46 should be used until the cells reach 80–90% confluence.

47

48 **▲ Note;**

- 49 - Cells should be passed until the cells reach 80–90% confluence.
- 50 - Cells should be used within 2 passages.

1
2 *Preparation of test chemicals, choice of solvent and precipitation evaluation*
3

4 Test chemicals are dissolved in an appropriate solvent such as PBS(-) and
5 DMSO. The recommended maximum final concentrations are 1 % (v/v) for
6 PBS(-) or 0.1 % (v/v) for DMSO. The chemicals must be weighed and
7 dissolved in solvent before each experiment, including 5-FU for the
8 positive-reference control. The highest test concentration of any chemicals is
9 1000 µg/mL.

10 Dilution series of 7 concentrations of 5-FU (0.0003, 0.0006, 0.0013, 0.0025,
11 0.005, 0.01 and 0.02 mg/mL; with a common ratio of 1:2) as positive-reference
12 control is dissolved from a 0.02 mg/mL prepared stock solution in PBS(-).
13

14 A flowchart for selection of appropriate vehicle, preparation of chemicals and
15 precipitation evaluation is shown in Figure 2.

16 Dissolve the test chemical first in PBS(-) to obtain a concentration of 100
17 mg/mL. For example, weigh 10 mg of the test chemical in an appropriate tube
18 and add PBS(-) up to 100 µL. If the chemical is not soluble at 100 mg/mL,
19 dilute with a common ratio of 1:2 to obtain a concentration of 50 mg/mL. If the
20 chemical is not soluble at 50 mg/mL, dilute with a common ratio of 1:2 to
21 obtain a concentration of 25 mg/mL. Other final volumes are also possible but
22 the minimum required one should be 100 µL depending on the needs of test
23 chemical.

24 If the chemical is not soluble at 25 mg/mL in PBS(-), the chemical should be
25 dissolved in DMSO at 1000 mg/mL. For example, weigh 100 mg of the test
26 chemical in an appropriate tube and add DMSO up to 100µL. At this step, a
27 precipitation test should be undertaken as followed: add 1µl of the
28 DMSO-dissolved chemical in the tube containing 999µl of assay medium. Mix
29 the medium with the chemical thoroughly by vortex. If the precipitate is not
30 observed at 1000µg/ml, ~~check water solubility over 1000 µg/ml.~~ then 1000
31 µg/ml should be the highest concentration tested.
32

33 If the chemical is not soluble at 1000 mg/mL, the highest soluble
34 concentration should be determined by diluting the solution from 500 mg/ml at
35 a common ratio of two (250 mg/ml -> 125 mg/ml if needed) with DMSO. If the
36 chemical is not soluble at 250 mg/mL in DMSO, prepare solution of test
37 chemical using PBS(-) or DMSO with 2-fold dilution factor (ex; 125 mg/mL in
38 DMSO or 12.5 mg/mL in PBS(-) if needed). Other final volumes are also
39 possible but the minimum required one should be 100 µL depending on the
40 needs of test chemical.


41 Sonication and vortexing may be used if needed, and attempt to dissolve the
42 chemical for at least 5 minutes.

43 For each highest concentration where the chemical is soluble, the
44 precipitation test should be provided. For example, if the chemical is soluble in
45 DMSO at 500 mg/mL, then from this concentration and with a 2-fold dilution
46 factor, the possibility of precipitation in the assay medium should be verified.
47 For this purpose, add 1 µL of the DMSO dissolved chemical in the tube
48 containing 999 µl of assay medium. Mix the medium with the chemical
49 thoroughly by vortex.

50 If a precipitate ~~appears on the tips of the 8 channel micropipette or~~ remains in

1 the medium after mixing, then the corresponding concentrations should not be
 2 taken into account for the experiment. For example, if the last precipitate is
 3 observed at 250 μ g/mL then the highest concentration used for the experiment
 4 will be 125 μ g/mL.

5 For chemicals that do not dissolve in PBS, the solubility in DMSO should be
 6 carefully evaluated especially for liquid chemicals where small and difficult to
 7 see bubbles may appear.
 8



9
10
11 Figure 2 Vehicle selection and preparation of test chemicals
12
13

14 *Serial dilution of test chemicals and arrangement of plate format*

15 Generally, the test chemicals should be prepared as follows;

- 16 - Maximum concentration (final concentration): 1000 μ g/mL (depend on the
- 17 solubility of test chemicals)
- 18 - Serially dilute with the same solvent at a common ratio of 1:10, basically in
- 19 the dose finding assay. Carry out the definitive study with smaller dilution
- 20 factor (ex; 1:1.5, 1:2, 1:3 or 1:5) covering the relevant range of dose
- 21 response and as much as possible and if the gap between the IC₅₀ and
- 22 ID₅₀ allows it, the lowest common ratio should be chosen. If IC₅₀ and ID₅₀
- 23 values are close to each other, the IC₅₀ value should be situated in the
- 24 middle of the new concentration range.
- 25
- 26

- 27 1. Prepare dissolved 5-FU as positive-reference control and test
- 28 chemicals according to the procedure (described in section
- 29 "Preparation of test chemicals") using 96 well microplates (Corning,
- 30 Cat. 3359) (cf. Figure 3, row 1A - 1H, 2A - 2H) before assay.

1

	Positive control (5-FU)		Test chemicals										
	1	2	3	4	5	6	7	8	9	10	11	12	
A	PBS(-)	VC											
B	0.0003 mg/ml	Conc. 1											
C	0.0006 mg/ml	Conc. 2											
D	0.0013 mg/ml	Conc. 3											
E	0.0025 mg/ml	Conc. 4											
F	0.005 mg/ml	Conc. 5											
G	0.01 mg/ml	Conc. 6											
H	0.02 mg/ml	Conc. 7											

VC: Vehicle control (DMSO or PBS(-))

Figure 3 Serial dilution and arrangement of test chemicals (1)

- Add 20 μ L of diluted 5-FU to the well of 96 well Assay Block containing 980 μ L of assay medium using an 8 channel micropipette (cf. Figure 4, row 1A – 1H).
- Add 20 μ L of chemicals to the well of 96 well Assay Block containing 980 μ L of assay medium using an 8 channel micropipette if chemicals are dissolved in PBS(-) (cf. Figure 4, row 2A – 2H).

Add 2 μ L of chemicals to the well of 96 well Assay Block containing 998 μ L of assay medium using an 8 channel micropipette if chemicals are dissolved in DMSO. Store solution at room temperature (20 - 30°C) until next steps.

	1	2	3	4	5	6	7	8	9	10	11	12	
A													
B													
C													
D													
E													
F													
G													
H													


Figure 4 Serial dilution and arrangement of test chemicals (2)

Dose finding study

Preparation procedure of cells for the assays.

- Dispense 100 μ L of PBS(-) into the peripheral wells and 100 μ L of assay medium into row 11B – 11G (cf. Figure 5).
- Prewarm assay medium and PBS(-) at 37 °C. 0.25% Trypsin /1mM EDTA solution is kept at room temperature (20 - 30°C) by bathing in water or ambient air.

- 1 3. Remove the supernatant from culture dish, wash with 5 mL of PBS(-).
 2 4. Add 1 - 2 mL of 0.25% trypsin/1 mM EDTA and remove immediately.
 3 Then, cells are incubated at 37 °C for 1 - 2 min.
 4 5. Add to 2 mL of assay medium to the dish, and suspend the cells using
 5 a micropipette (1000 μ L).
 6 6. Use cell strainer with 50mL centrifuge tubes to obtain more uniform
 7 single-cell suspensions.
 8 7. Check viability by staining an aliquot of the cell suspension with
 9 Trypan blue.
 10 **▲Note** A viability of 90 % or higher is acceptable.
 11 8. Count cell number.
 12 9. Dilute the cells in assay medium at cell density of 15,000 cells/mL (750
 13 cells /50 μ l/well).
 14 10. Transfer the cells to a reservoir, and dispense 50 μ L of cell suspension
 15 to PrimeSurface® 96well white plate (Sumitomo Bake Lite, Co. Ltd.,
 16 MS-9096W) (cf. Figure 5; row 2B – 10G in the plates for 5-FU or test
 17 chemicals).
 18 11. Incubate the cells at 37°C with 5% CO₂.
 19
 20



- 21
 22
 23
 24
 25
 26
 27
 28
 29
 30
 31
 32
 33
 34
 35
 36
 37

Figure 5 Preparation of the cells for the assay

Treatment of the test chemicals for the assays

1. After more than 2 hrs of incubation, mix prepared solution of 5-FU and test chemicals by pipetting thoroughly.
2. Add 50 μ L of assay medium to a plate (cf. Figure 6; row 2B –2G).
3. Dispense 50 μ L of solution of chemicals to 96 well-plate containing the cells using an 8 channel micropipette (cf. Figure 6; row 3B –10G).
4. Shake the plate with a plate shaker for a few seconds.
5. Incubate the cells at 37°C with 5% CO₂.




Figure 6 Arrangement of chemicals and vehicles for assay

Measurement of fluorescence and luminescence

▲ Note Switch on measuring equipments 30 min before starting the measurement.

120 hours after seeding the Hand1-ES cells, measure the two following endpoints as described below:

For cytotoxicity assay

1. Completely thaw the CellTiter-Fluor™ Cell Viability Assay components in a 37°C water bath. Vortex the GF-AFC substrate (**▲ Note** avoid exposure to light) to ensure homogeneity, then briefly centrifuge for complete substrate volume recovery.
2. Transfer the GF-AFC Substrate into the Assay Buffer container to form a complete reagent (cf; Assay buffer 1mL : GF-AFC substrate 10 μL) . Mix by vortexing the contents until the substrate is thoroughly dissolved.
3. Add 10 μL of CellTiter-Fluor™ complete reagent to each well (cf. Figure 6; row 2B - 11G).
4. Mix briefly by orbital shaking, then incubate for at least 30 minutes at 37°C. Seal the surface of the plates.
▲ Note Don't incubate longer than 1hour, and be sure to shield plates from ambient light.
5. Measure resulting fluorescence using a fluorometer (380–400nmEx /505nmEm).
▲ Note You may need to adjust instrument gains.

For differentiation assay

1. Thaw the Steady-Glo® Luciferase Assay System and equilibrate to 25 °C prior to use (**▲ Note** avoid exposure to light).
2. Dispense 100 μL of Steady-Glo® Luciferase Assay System to each well of 96 well plates by an 8 channel micropipette (cf: Figure 6; row 2B – 11G).
3. Shake the plates gently for over 30 min at room temperature (20 – 30 °C) with a plate shaker. Seal the surface of the plates.

- 1 **▲ Note** Don't incubate longer than 1hour, and be sure to shield plates
 2 from ambient light.
 3 4. Place the plate in luminometer. Luminescence is measured with the
 4 appropriate time (cf. 1sec for validation study) according to the
 5 manufacturer's instructions.

6
7
8 *Definitive study*

9 Calculate endpoints according to Steps "Data analysis". For the definitive
 10 study, choose seven dilutions. Carry out the definitive study with a smaller
 11 dilution factor covering the relevant range of dose response according to
 12 dose finding study ([see Section of Experimental validity of definitive](#)
 13 [studies](#))

14
15
16 *Data analysis*

17 Use an Excel (Microsoft) spreadsheet for data recording (Result format ver.
 18 08E). Calculation method and the acceptance criteria for these assays are
 19 as follows;

20
21 Determination of ID₅₀ and IC₅₀ values

22
23 Relative viabilities (surviving rate) or relative activities for each
 24 concentration of test chemicals against vehicle control can be calculated
 25 using Excel (Microsoft) spreadsheet. Calculation formulas are described as
 26 follows;

27
28 For differentiation assay

- 29 a : Mean values of vehicle control
 30 b : Mean values of test chemicals
 31 c : Mean values of background
 32 d : NET values of vehicle control : a-c
 33 e : NET values of test chemicals : b-c
 34 A : Relative activity of vehicle control : 100 (%)
 35 B : Relative activity of test chemicals : e/d X 100 (%)

36
37 For cytotoxicity assay

- 38 a : Mean values of vehicle control
 39 b : Mean values of test chemicals
 40 c : Mean values of background
 41 d : NET values of vehicle control : a-c
 42 e : NET values of test chemicals : b-c
 43 A : Relative viability of cells for vehicle control : 100 (%)
 44 B : Relative viability of cells for test chemicals : e/d X 100 (%)

45
46 Inhibition of differentiation is expressed as the concentration of the test
 47 chemical that reduced the luminescence by 50% (ES-ID₅₀, calculated from
 48 the concentration-response curve). Cytotoxicity is expressed as the
 49 concentration of the chemical reducing the viability of cells to 50% of the
 50 control level (ES-IC₅₀ determined from concentration-response curves).

51 The IC₅₀ and ID₅₀ values relative to solvent controls were derived from

1 computational regression analysis of concentration response curves with a
2 three parameter logistic function using an Excel (Microsoft) spreadsheet for
3 data recording (Result format ver.08E).

4 The full concentration-response curve is required for the calculation of the
5 ID_{50} and IC_{50} values, but this may not always be achievable or practical due
6 to limitations of the test concentration range (for example due to cytotoxicity
7 or solubility problems).

8
9
10 *Acceptance criteria*

11
12 **Quality control**

13 **i) Quality control of cell condition**

14 Check viability of the cells by staining an aliquot of the cell suspension with
15 Trypan blue. A viability of 90% or greater is acceptable.

16
17 **ii) Quality control of differentiation and cytotoxicity assays**

18 To verify the cell growth and cell differentiation, the comparison of the
19 medium control (MC) and background (BG) is used.

20 - The Lower limit of the 95% confidence intervals of the ratio of MC / BG
21 should be above 1 for the cytotoxicity assay

22 - The Lower limit of the 95% confidence intervals of the ratio of MC / BG
23 should be above 10 for differentiation assay

24
25 **iii) Performance standard of the assay**

26 The quality of the assay must be controlled using 5-FU as a
27 positive-reference chemical.

28 -The range of ID_{50} for 5-FU should be within **0.003 and 0.067 $\mu\text{g/mL}$**
29 (according to phase 0 study).

30 -The range of IC_{50} for 5-FU should be within **0.003 and 0.065 $\mu\text{g/mL}$**
31 (according to phase 0 study).

32
33 **iv) Quality control for effect of vehicle**

34 To verify the effect of the vehicle, the comparison of the medium control
35 (MC) and vehicle control (VC) is used.

36 - The Lower limit of the 95% confidence intervals of the ratio of VC / MC
37 should be above 0.2 for the cytotoxicity and differentiation assays.

38 - In addition to the above criterion

39 CV (SD/mean X 100%) of VC should be below 100%.

~~1 EXPERIMENT 3 (for optional prediction with PM2)~~

~~2 3-1 Measurement of metabolic stability of test chemicals~~

~~3 4 *Liver S9 fractions*~~

~~5 Rat liver S9 fractions are commercially available (Male, pool of 400 livers, XenoTech). The brief procedure for preparation is as follows. A portion of liver from male SD rats is homogenized in 4 volumes of Tris/HCl buffer (pH7.4) with a homogenizer. Liver homogenate is centrifuged at 9,000 × g for about 20 min at 4 °C to separate S9 fractions. The protein level in the S9 fraction can be measured by commercially available protein assay kit.~~

~~6 12 *Metabolic stability*~~

~~7 Substrate solution is prepared by mixing 100 μM substrate in acetonitrile with 100 fold volume of 6 mM β-NADPH (Oriental Yeast) in 125 mM potassium phosphate buffer (pH7.4). The reaction is initiated by mixing 50 μL of substrate solution with 50 μL of diluted liver S9 fractions (0.5 mg protein/mL) in 125 mM potassium phosphate buffer (pH7.4). After incubation for 35 min at 37 °C, the reaction was stopped by addition of 400 μL methanol. After centrifugation at 3,000 g for 10 min, the supernatant is analyzed by the triple quadrupole LC/MS and the peak area of the target ion is calculated. The control sample is prepared without liver S9 fractions. The metabolic stability is expressed as the ratio of peak area with liver S9 to that of control.~~

JUDGMENT

In the Hand1-Luc EST, the embryotoxic potency of test chemicals was predicted using following equation. ~~a biostatistically based prediction model shown in Figure 7.~~

Discriminant equation

Prediction Model 1(PM1)

$$\text{Score} = 1.474 \times \log(1) + 0.714 \times \log(2) - 1.109$$

$$1, \log(\text{IC}_{50}/\text{ID}_{50}) ; 2, \log(\text{Maximum dose}/\text{IC}_{50})$$

Probability: inverse logit transform of the score

Positive : probability ≥ 0.52

Negative : probability < 0.52

Majority judgment of 3 replicates (~~Phase 2a~~)

	Dose finding	Definitive study1	Definitive study2	Definitive study3	Final judgment
Case 1	-	P	P	-	P
Case 2	-	N	P	P	P
Case 3	-	N	P	N	N
Case 4	-	N	N	-	N
Case 5	-	P	N	P	P
Case 6	-	P	N	N	N

The dose finding is not considered for the evaluation. If both first definitive studies conclude negativity or positivity, then only two definitive studies are required.

However, if definitive study 1 and definitive study 2 have 2 different conclusions, then a 3rd definitive study will be required.

Experimental validity of definitive studies:

If the IC₅₀ and ID₅₀ values of DS1 compared to the IC₅₀ and ID₅₀ values of DS2 are different from a factor 7, then one more test should be done and the 2 closest results should be taken into account.

(Example1: If DS1, IC₅₀=10 µg/ml, DS2, IC₅₀=125 µg/ml and DS3, IC₅₀=28 µg/ml then only the results of DS1 and DS3 should be analyzed.

Example 2: If DS1, IC₅₀=10 µg/ml, DS2, IC₅₀=12 µg/ml and DS3, IC₅₀=128 µg/ml then only the results of DS1 and DS2 should be analyzed)

If, again, IC₅₀ and ID₅₀ values of DS3 are higher or lower than a factor 7 compared to DS1 and DS2, then one more last test should be performed. This evaluation should be done disregarding the positivity or negativity judgment.

1
2
3
4
5

OPTION (Experiment 2)

In the Hand1-Luc EST, prediction of the embryotoxic potency of test chemicals can be improved optionally by measuring metabolic stability (to be submitted).

UPDATE RECORD

- 1
2
3
4 Ver.01E, 2013, Feb, 20th distributions
5
6 Ver.02E, 2013, April, 30th distributions
7 Major modification
8 - Plate format for the assays
9 - Positive reference chemicals
10 - Preparation of test chemicals
11 - Maximum final concentrations of vehicle (PBS(-) and DMSO)
12 - Quality control and acceptability criteria
13
14 Ver.03E, 2013, July, 18th distributions
15 - Plate format for the assays
16 - Detection method for cytotoxicity
17
18 Ver.04E, 2013, July 23th distributions
19 - Minor modification
20 - Tentative acceptance criteria
21
22 Ver.05-1 E, 2013, November 6th distributions
23 Major modification
24 - Measurement at 120 hr
25 - Tentative acceptance criteria or phase 1 study.
26
27 Ver.06-1 E, 2014, March 10th distributions
28 Major modification:
29 - Measurement of the water solubility, suppression of the measurement
30 of the logP
31 - Inclusion of a precipitation evaluation step
32 - Modification of the prediction model
33 - Tentative acceptance criteria modified for phase 2.
34
35 Ver.07 E, 2014, May 14th distributions
36 Minor modification:
37 - Measurement of the water solubility
38 - Modification of the prediction model
39
40 Ver.08 E, 2015, June 26th distributions
41 Minor modification:
42 - According the suggestion of VMT member, calculation method of IC₅₀
43 and ID₅₀ values using a three parameter logistic function instead of
44 2-parameter
45 - Modification of the prediction model
46
47 Ver.09 E, 2016, March 9th distribution
48 Minor modifications:
49 - Correction English figure 1
50 - Correction on the top seal product reference number (expendable
51 supplies section), that has been updated by the supplier.

- 1 - Precisions about dissolution of liquid chemicals in solvent
- 2 - Figure 2 corrected with insertion of evaluation of precipitation in the
- 3 assay medium
- 4 - Precision on dilution ratio that should be chosen for the definitive
- 5 studies
- 6 - Definitive studies' quality control: IC₅₀ and ID₅₀ of DS1 and the
- 7 ones in DS2 are significantly different.
- 8 - Correction of the prediction model.
- 9 - Correction of optional experiment for PM2