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1）OECD (2025), Test Guideline No. 442E *In Vitro* Skin Sensitisation assays addressing the Key Event on activation of dendritic cells on the Adverse Outcome Pathway for Skin Sensitisation

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| | |
|---|---|
| DAPI: | 4',6-Diamidino-2-phenylindole |
| GAPDH: | Glyceraldehyde-3-phosphate dehydrogenase |
| GHS: | Globally Harmonized System |
| h-CLAT: | Human Cell Line Activation Test |
| IL-8 Luc assay: | Interleukin-8 Reporter Gene Assay |
| LLNA: | Local Lymph Node Assay |
| OECD: | Organisation for Economic Co-operation and Development |
| PI: | Propidium iodide |
| TG: | Test Guideline |
| U-SENS™: | U937 Cell Line Activation Test |

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| パラメーター | 定義 |
|---|---|
| GAPLA | SLR luciferase activity driven by GAPDH promoter |
| IL8LA | SLO luciferase activity driven by IL-8 promoter |
| nIL8LA | IL8LA/GAPLA |
| Ind-IL8LA | nIL8LA of THP-G8 cells treated with chemicals / nIL8LA of untreated cells |
| Inh-GAPLA | GAPLA of THP-G8 cells treated with chemicals / GAPLA of untreated cells |
| PI reduction | % of PI-excluding cells among THP-G8 cells treated with chemicals / % of PI-excluding cells among THP-G8 cells without chemical treatment. |

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１）OECD (2025), Test Guideline No. 442E *In Vitro* Skin Sensitisation assays addressing the Key Event on activation of dendritic cells on the Adverse Outcome Pathway for Skin Sensitisation

２）OECD (2023), Validation report for the international validation study on the IL-8 Luc assay as a test evaluating the skin sensitizing potential of chemicals conducted by the IL-8 Luc Assay. Series on Testing and Assessment No. 267. Organisation for Economic Cooperation and Development, Paris. Available at: (https://one.oecd.org/document/ENV/CBC/MONO(2017)19/REV1/en/pdf)

３）Kimura, Y., Fujimura C. and Aiba, S. (2021) The modified IL-8 Luc assay, an in vitro skin sensitisation test, can significantly improve the false-negative judgment of lipophilic sensitizers with logKow values > 3.5 Archives of Toxicology 95(3):1-10, DOI:10.1007/s00204-020-02934-9

４）OECD (2025) OECD Guideline No. 497. Defined Approaches on Skin Sentitisation. Organization for Economic Cooperation and Development, Paris. Available at: (oecd-ilibrary.org)